




SKU: S0245

SHOWER COMMODE - TILT IN SPACE MOBILE - 445 MM

From: \$3,999.99

The Tilt in Space Mobile Shower Commode assists in the transportation and care for patient showering and toileting. This commode comes with extra high PVC strapped backrest for improved comfort and safety. Also comes with an integrated push-handle for manoeuvrability and has a 10 -35 degree tilt angle. Equipped with swing-away footrest and contoured padded calf supports. This commode has been independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind .

Please note, seat sold separately



Scan QR code to
view product online

PRODUCT DESCRIPTION

- Independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind
- Extra high PVC strap backrest for improved comfort and safety
- 10-35° tilt angle
- Integrated push-handle for manoeuvrability
- Swing-up padded arms for unrestricted side transfer
- Swing-away and removable height adjustable footrest with contoured padded calf supports
- 4x 125mm total locking castors
- Available in 400, 445, 500, 550 and 600mm seat widths
- Padded seat not included

ADDITIONAL INFORMATION

Overall Width	580 mm
Overall Depth	1070 mm
Overall Height	1453 mm
Seat Width	445 mm
Seat Depth	430 mm
Floor To Seat Height	540 mm
Armrest Type	Padded Loop
Armrest Height	365 mm
Backrest Height	855 mm
Footrest Type	Swing Away
Product Weight	21 kg
Max User Weight	200 kg
Brand	KCare