



SKU: S0245

SHOWER COMMODE - MOBILE - TILTING

From: \$3,999.99

The Tilt in Space Mobile Shower Commode assists in the transportation and care for patient showering and toileting. This commode comes with extra high PVC strap backrest for improved comfort and safety. An integrated push-handle for manoeuvrability and has a 10 -35 degree tilt angle. Equipped with swing-away footrest and contoured padded calf supports.

Seat sold separately



Scan QR code to
view product online

PRODUCT DESCRIPTION

Key Features

- Independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind
- Extra high PVC strap backrest for improved comfort and safety
- 10-35° tilt angle
- Integrated push-handle for manoeuvrability
- Swing-up padded arms for unrestricted side transfer
- Swing-away and removable height adjustable footrest with contoured padded calf supports
- 4x 125mm total locking castors
- Designed to be used with a K Care patented removable seats
- Available in 400, 445mm and 500mm seat widths
- Padded seat not included.

ADDITIONAL INFORMATION

Brand	KCare
Overall Width	580 mm
Overall Depth	1070 mm
Overall Height	1453 mm
Seat Width	445 mm
Seat Depth	430 mm
Floor To Seat Height	540 mm
Armrest Type	Loop
Arm Height	365 mm
Backrest Height	855 mm
Footrest Type	Swing Away
Product Weight	21 kg
Max User Weight	200 kg