



**SKU:** S0245

## **SHOWER COMMODE - TILT IN SPACE MOBILE - 445 MM**

**From: \$3,999.99**

The Tilt in Space Mobile Shower Commode assists in the transportation and care for patient showering and toileting. This commode comes with extra high PVC strapped backrest for improved comfort and safety. Also comes with an integrated push-handle for manoeuvrability and has a 10 -35 degree tilt angle. Equipped with swing-away footrest and contoured padded calf supports. This commode has been independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind .

**Please note, seat sold separately**



Scan QR code to  
view product online

## **PRODUCT DESCRIPTION**

- Independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind
- Extra high PVC strap backrest for improved comfort and safety
- 10-35° tilt angle
- Integrated push-handle for manoeuvrability
- Swing-up padded arms for unrestricted side transfer
- Swing-away and removable height adjustable footrest with contoured padded calf supports
- 4x 125mm total locking castors
- Available in 400, 445, 500, 550 and 600mm seat widths
- Padded seat not included

## ADDITIONAL INFORMATION

<b>Overall Width</b>	580 mm
<b>Overall Depth</b>	1070 mm
<b>Overall Height</b>	1453 mm
<b>Seat Width</b>	445 mm
<b>Seat Depth</b>	430 mm
<b>Floor To Seat Height</b>	540 mm
<b>Armrest Type</b>	Padded Loop
<b>Armrest Height</b>	365 mm
<b>Backrest Height</b>	855 mm
<b>Footrest Type</b>	Swing Away
<b>Product Weight</b>	21 kg
<b>Max User Weight</b>	200 kg
<b>Brand</b>	KCare