



SKU: S7914

DRIVE MEDICAL TITAN POWER CHAIR

\$3,299.99

The Drive Medical Titan Powered Wheelchair is a stylish, front-wheel drive mobility solution suited to both indoor and outdoor spaces. With a 63.5 cm turning radius and tool-free disassembly, it offers flexibility, transportability and comfort.



Scan QR code to
view product online

GALLERY IMAGES





PRODUCT DESCRIPTION

The Drive Medical Titan Powered Wheelchair combines portability with versatile performance for those who need a reliable indoor and outdoor power wheelchair.

Its front-wheel drive configuration improves stability and traction, enabling the chair to manage small obstacles and uneven terrain with ease.

With a 63.5 cm turning circle, the Titan navigates confined spaces, making it suitable for use at home or in busy public settings.

The chair disassembles without tools, allowing loading into vehicles for transport. Its compact design includes a folding backrest and a swivel seat for easy transfers.

Comfort is supported through padded, height-adjustable armrests, an adjustable and reclining seat, and a removable headrest.

A dynamic Shark controller provides smooth handling, while solid, non-marking tyres ensure durability and minimal upkeep.

Safety features include a seat belt and de-clutch functionality for controlled freewheeling when required.

Who it's for

- For clients who need a powered wheelchair that is easy to transport and suitable for varied environments.
- For therapists seeking a stable, manoeuvrable option with multiple comfort adjustments for their client.

Key features

- **Front-wheel drive system** – enhances stability and outdoor performance
- **63.5 cm turning radius** – supports use in tight indoor spaces

- **Tool-free disassembly** – for simple transport and storage
- **Folding backrest and swivel seat** – assist with transfers and compact storage
 - **Height-adjustable, reclining seat and headrest** – for custom comfort
 - **Detachable, padded armrests** – adapt to individual needs
- **Dynamic Shark programmable controller** – provides smooth and intuitive driving
 - **Solid, non-marking tyres** – reduce maintenance and improve durability
 - **De-clutch function** – allows safe freewheeling
 - **Seat belt included** – enhances safety during use

ADDITIONAL INFORMATION

Safe Working Load (SWL)	136 kg
Battery Type	Sealed Lead Acid
Battery Size	2 x 12V/35 AH
Range	24 km
Max speed	6.5 km/h
Size (W x L assembled)	63.5 cm x 102 cm
Seat Width	45.7 cm
Seat Depth	45.7 cm
Seat Height	55 cm
Weight	72 kg
Front Wheels	10" Drive Wheels
Rear Wheels	6" Castor Wheels
Max Safe Incline Slope	6°
Turning Circle	63.5 cm
Brand	Drive Medical

[Drive Medical Titan Brochure](#)

[Drive Medical Titan User Manual](#)

Talk to us on the phone

Our experienced phone sales team can assist you with any product enquiries and help you make a purchase directly on the phone. We are available anytime between 9am and 5pm, Monday to Friday.

[1300 008 267](tel:1300008267)

ARE THE TYRES INFLATABLE OR SOLID RUBBER?

The tyres are **solid puncture proof** and it has anti-tip wheels at the front as standard.