




SKU: S7914

DRIVE MEDICAL TITAN POWER CHAIR

\$2,999.99

The Titan is also great indoors, navigating tight, enclosed spaces with ease.

- Designed for indoor and outdoor use
Front wheel drive system for optimum manoeuvrability
- Detachable padded and adjustable armrests to suit the comfort of the user
- Adjustable length controller mount
- Folding backrest for easy storage and transportation
- Height adjustable and reclining seat to suit comfort preferences
- Swivel seat facilitates easy transfers on and off the chair
- Height adjustable and removable headrest to suit user requirements
- The de-clutch operation to enable you to safely freewheel the chair where and when required
- Padded upholstery for additional comfort and support



Scan QR code to
view product online

PRODUCT DESCRIPTION

The front wheel drive Titan Power chair is a comfortable and stylish option for users seeking a transportable power mobility product that's ideal for indoor and outdoor use. The Titan's front-wheel system positions the drive wheel forward of the user's centre of gravity, making the chair extremely stable and giving it great performance over small obstacles.

ADDITIONAL INFORMATION

Brand	Drive Medical
Width	63.5cm
Length	102cm
Height	112cm
Seat Width	46cm
Seat Depth	46cm
Seat Height	55cm
Weight	72kg
Max SWL (Maximum Weight)	136kg
Heaviest Piece Disassembled	25 kgs
Seat Weight	19.35 kgs
Battery Weight (PAIR)	19.9 kgs.
Base Weight	25kg's
Charger	3A Offboard
Maximum speed	6.48 km/h
Maximum range	24 km
Width between the armrests	43cm - 53cm
Warranty	2 years