



**SKU:**

## **ROYALE TODD ATLANTIC HIGHBACK DAY CHAIR**




**\$1,049.99**

The Royale Todd Atlantic High Back Day Chair provides robust, adjustable bariatric seating. Designed for hospital, rehabilitation and home environments, it supports posture, comfort and daily use.



Scan QR code to  
view product online

## VARIATIONS

Image	SKU	Price	Description	Colour
	S13876	\$1,049.99		Beige Vinyl
	S13878	\$1,049.99		Light Grey Fabric
	S13879	\$1,049.99		Storm Grey Vinyl

## PRODUCT DESCRIPTION

The Royale Todd Atlantic High Back Day Chair is a versatile and supportive bariatric seating option for those needing enhanced comfort.

It features adjustable seat height and depth, allowing a customised fit for different users and clinical needs.

The contoured backrest offers ergonomic lumbar support, while the adjustable neck pillow provides targeted comfort for the head and neck. The padded seat, armrests and backrest enhance comfort during prolonged sitting.

Upholstered in durable fabric or vinyl, the chair is easy to clean and suited for both home and care settings.

The powder-coated steel frame adds strength and corrosion resistance, supporting individuals up to 350 kg.

### Who it's for

- For clients who need high back bariatric seating with customisable support and comfort.
- For therapists seeking a high-capacity, adjustable day chair for bariatric clients across clinical and home settings.

### Key features

- **Adjustable seat height and depth** – enables a personalised fit for different body types
- **Contoured high backrest with lumbar support** – promotes upright posture and comfort
- **Adjustable neck pillow** – supports the head and neck for added relaxation
- **Vinyl or fabric upholstery options** – suit varied environments and hygiene needs
- **Powder-coated steel frame** – enhances durability and resists corrosion
- **High weight capacity** – safe working load of up to 350 kg

## ADDITIONAL INFORMATION

<b>Brand</b>	Royale Medical
<b>Seat Width</b>	500 mm
<b>Seat Depth</b>	420 – 500 mm
<b>Seat Height</b>	420 – 600 mm
<b>Arm Height</b>	220 mm
<b>Back Height</b>	740 mm
<b>Overall Width</b>	670 mm
<b>Max SWL</b>	350 kg
<b>Colour</b>	Beige Vinyl, Cream Fabric, Light Grey Fabric, Storm Grey Vinyl