




SKU:

TODD ATLANTIC CHAIR

\$1,049.99

The Todd Atlantic Adjustable Lounge chair is now one of the most widely accepted chairs for Rehabilitation or Hospital seating, where the ability to adjust the shape of the chair to suit a variety of patients and their needs is important. The shaped seat and backrest, well positioned armrests and ease of adjustment make Todd Atlantic the chair of choice for leading Hospitals and Rehab organisations.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Colour
	S13876	\$1,049.99		Beige Vinyl
	S13878	\$1,049.99		Storm Grey Fabric
	S13879	\$1,049.99		Storm Grey Vinyl

PRODUCT DESCRIPTION

The seat length is easily adjusted by loosening the knobs beneath the seat, sliding the backrest out to create the required seat length and re-tightening the knobs. The seat height is adjusted by pressing the adjuster pin and moving the leg to the desired length, and the adjustable neckrest can be positioned to suit.

To increase comfort and enhance posture, the backrest has built-in lumbar, lateral and head support. The seat is shaped with a front bolster to resist extension and assist in posture enhancement. The adjustable neck rest Dacron padding can be altered or removed.

Features

- High back orthopaedic chair
- Adjustable height and seat depth
- Contoured buttoned backrest provides added lumbar and neck support
- Adjustable neck pillow
- Concealed back frame for hygiene and infection control purposes
- Powder coated steel frame

ADDITIONAL INFORMATION

Brand	Royale Medical
Seat Width	500mm
Seat Depth	420-500mm
Seat Height	420-600mm
Arm Height	220mm
Back Height	740mm
Overall Width	670mm
Max SWL	350kg
Colour	Beige Vinyl, Cream Fabric, Storm Grey Fabric, Storm Grey Vinyl