

SKU: S2439




TOPFORM ASHLEY ELECTRIC RECLINER LIFT CHAIR MAXI

~~\$2,499.99~~ \$2,499.99

Australian made, the Topform Ashley Lift Chair is a luxury recliner chair designed with great features including a Tasmanian Oak frame with four optional (37mm) castors, high density foam padded seat and backrest with 3 separate padded horizontal cushions. Wide variety of colour/fabric choices available.

For colour choices please email
onlineorders@ilsau.com.au



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Motor
	S2443	\$2,499.99		Twin
	S2439	\$2,099.99		Single

PRODUCT DESCRIPTION

The Topform Ashley Lift Chair is a luxurious recliner chair with a Tasmanian Oak frame and four optional castors.

Full length, fully padded and upholstered armrests and available in a standard range of fabrics, non-standard fabrics at an additional cost. A two year frame, motor, transformer and hand control warranty and 12 month springs and component warranty are standard across the Topform Lift Chair range.

The lifting mechanism has been specifically designed to assist those who need help to get out of the chair with ease and without any body strain. Available as a single motor or twin motor lift chair.

Single vs twin vs quad motor lift chairs - what's the difference?

All ILS lift chairs have a lifting action, to support the user to stand from a seated position.

- Single motor chairs move the back rest and leg rest at the same time.
- Twin motor chairs allow the backrest and leg rest to move independently of each other, typically allowing for a greater degree of leg elevation and greater degree of recline compared to a standard single motor chair.
- Quad motor chairs allow independent movement of headrest, lumbar, backrest and footrest allowing for the greatest level of re-positioning of all chairs.

ADDITIONAL INFORMATION

Seat Height	530 mm
Seat Depth	520 mm
Seat Width	530 mm
Seat To Top Of Back	680 mm
Chair Width	830 mm
Chair Height	1170 mm
Max SWL	155 kg
User Height	180 cm to 188 cm
Brand	Topform
Note	Made in Australia is various fabrics, vinyls and leathers. All Topform chairs are custom made to order and involve production lead times*