

SKU: S2439




## TOPFORM ASHLEY ELECTRIC RECLINER LIFT CHAIR MAXI

**~~\$2,499.99~~ \$2,499.99**

Australian made, the Topform Ashley Lift Chair is a luxury recliner chair designed with great features including a Tasmanian Oak frame with four optional (37mm) castors, high density foam padded seat and backrest with 3 separate padded horizontal cushions. Wide variety of colour/fabric choices available.

For colour choices please email  
[onlineorders@ilsau.com.au](mailto:onlineorders@ilsau.com.au)



Scan QR code to  
view product online

VARIATIONS

| Image                                                                            | SKU   | Price      | Description | Motor  |
|----------------------------------------------------------------------------------|-------|------------|-------------|--------|
|  | S2443 | \$2,499.99 |             | Twin   |
|  | S2439 | \$2,099.99 |             | Single |

## PRODUCT DESCRIPTION

The Topform Ashley Lift Chair is a luxurious recliner chair with a Tasmanian Oak frame and four optional castors.

Full length, fully padded and upholstered armrests and available in a standard range of fabrics, non-standard fabrics at an additional cost. A two year frame, motor, transformer and hand control warranty and 12 month springs and component warranty are standard across the Topform Lift Chair range.

The lifting mechanism has been specifically designed to assist those who need help to get out of the chair with ease and without any body strain. Available as a single motor or twin motor lift chair.

### **Single vs twin vs quad motor lift chairs - what's the difference?**

All ILS lift chairs have a lifting action, to support the user to stand from a seated position.

- Single motor chairs move the back rest and leg rest at the same time.
- Twin motor chairs allow the backrest and leg rest to move independently of each other, typically allowing for a greater degree of leg elevation and greater degree of recline compared to a standard single motor chair.
- Quad motor chairs allow independent movement of headrest, lumbar, backrest and footrest allowing for the greatest level of re-positioning of all chairs.

## ADDITIONAL INFORMATION

|                            |                                                                                                                                                 |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Seat Height</b>         | 530 mm                                                                                                                                          |
| <b>Seat Depth</b>          | 520 mm                                                                                                                                          |
| <b>Seat Width</b>          | 530 mm                                                                                                                                          |
| <b>Seat To Top Of Back</b> | 680 mm                                                                                                                                          |
| <b>Chair Width</b>         | 830 mm                                                                                                                                          |
| <b>Chair Height</b>        | 1170 mm                                                                                                                                         |
| <b>Max SWL</b>             | 155 kg                                                                                                                                          |
| <b>User Height</b>         | 180 cm to 188 cm                                                                                                                                |
| <b>Brand</b>               | Topform                                                                                                                                         |
| <b>Note</b>                | Made in Australia is various fabrics, vinyls and leathers.<br>All Topform chairs are custom made to order and involve<br>production lead times* |