

SKU: S2437



## TOPFORM ASHLEY ELECTRIC RECLINER LIFT CHAIR MINI

~~\$2,499.99~~ **\$2,499.99**

The Topform Ashley Lift Chair is a luxurious recliner chair with a Tasmanian Oak frame and four optional (37mm) castors, as well as a high density foam padded seat and backrest with three distinct padded horizontal cushions. Made in Australia.

Please email [onlineorders@ilsau.com.au](mailto:onlineorders@ilsau.com.au) with your color preferences.



Scan QR code to  
view product online

### VARIATIONS

Image	SKU	Price	Description	Motor
	S2437	\$2,099.99		Single
	S2441	\$2,499.99		Twin

## PRODUCT DESCRIPTION

Armrests are full length, completely padded and upholstered, and come in a standard range of fabrics, with non-standard fabrics available for an extra fee. The Topform Lift Chair comes with a two-year frame, motor, transformer, and hand control guarantee, as well as a 12-month spring and component warranty. Available as a single motor or twin motor lift chair.

### **Single vs twin vs quad motor lift chairs - what's the difference?**

All ILS lift chairs have a lifting action, to support the user to stand from a seated position.

- Single motor chairs move the back rest and leg rest at the same time.
- Twin motor chairs allow the backrest and leg rest to move independently of each other, typically allowing for a greater degree of leg elevation and greater degree of recline compared to a standard single motor chair.
- Quad motor chairs allow independent movement of headrest, lumbar, backrest and footrest allowing for the greatest level of re-positioning of all chairs.

## **ADDITIONAL INFORMATION**

<b>Seat Height</b>	470 mm
<b>Seat Depth</b>	490 mm
<b>Seat Width</b>	480 mm
<b>Seat To Top Of Back</b>	600 mm
<b>Chair Width</b>	780 mm
<b>Chair Height</b>	1020 mm
<b>Max SWL</b>	155 kg
<b>User Height</b>	153 cm to 160 cm
<b>Brand</b>	Topform
<b>Note</b>	Made in Australia is various fabrics, vinyls and leathers. All Topform chairs are custom made to order and involve production lead times*