

SKU: S2437


TOPFORM ASHLEY ELECTRIC RECLINER LIFT CHAIR MINI

~~\$2,499.99~~ \$2,499.99



The Topform Ashley Lift Chair is a luxurious recliner chair with a Tasmanian Oak frame and four optional (37mm) castors, as well as a high density foam padded seat and backrest with three distinct padded horizontal cushions. Made in Australia.

Please email onlineorders@ilsau.com.au with your color preferences.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Motor
	S2437	\$2,099.99		Single
	S2441	\$2,499.99		Twin

PRODUCT DESCRIPTION

Armrests are full length, completely padded and upholstered, and come in a standard range of fabrics, with non-standard fabrics available for an extra fee. The Topform Lift Chair comes with a two-year frame, motor, transformer, and hand control guarantee, as well as a 12-month spring and component warranty. Available as a single motor or twin motor lift chair.

Single vs twin vs quad motor lift chairs - what's the difference?

All ILS lift chairs have a lifting action, to support the user to stand from a seated position.

- Single motor chairs move the back rest and leg rest at the same time.
- Twin motor chairs allow the backrest and leg rest to move independently of each other, typically allowing for a greater degree of leg elevation and greater degree of recline compared to a standard single motor chair.
- Quad motor chairs allow independent movement of headrest, lumbar, backrest and footrest allowing for the greatest level of re-positioning of all chairs.

ADDITIONAL INFORMATION

Seat Height	470 mm
Seat Depth	490 mm
Seat Width	480 mm
Seat To Top Of Back	600 mm
Chair Width	780 mm
Chair Height	1020 mm
Max SWL	155 kg
User Height	153 cm to 160 cm
Brand	Topform
Note	Made in Australia is various fabrics, vinyls and leathers. All Topform chairs are custom made to order and involve production lead times*