



SKU:

TOPFORM LUXOR OPTIMA LIFT CHAIR

~~\$3,999.99~~ \$3,999.99

The Topform Luxor Optima Lift Chair offers pressure-reducing support for those with low risk of pressure injury. Made in Tasmania, it features multi-layered foam and multi-stretch fabric for enhanced comfort, and is available in multiple motor and size configurations.



Scan QR code to
view product online

VARIATIONS





Image	SKU	Price	Description	Motor	Sizes
	S8330	\$3,999.99		Twin	Maxi
	S8329	\$3,999.99		Twin	Medium
	S8328	\$3,999.99		Twin	Mini
	S8327	\$3,699.99		Single	Maxi

Image	SKU	Price	Description	Motor	Sizes
	S8326	\$3,699.99		Single	Mini

PRODUCT DESCRIPTION

The Topform Luxor Optima Lift Chair is built for comfort and low-risk pressure management. Designed and manufactured in Tasmania, it suits a wide range of individuals and is available in single or twin motor configurations, with three sizes to match individual needs.

Its multi-layered, pressure-reducing foam seat is covered in a breathable, two-way stretch material that is waterproof and anti-bacterial. This helps to promote airflow and reduce shear, supporting skin integrity. The Luxor Optima's lift mechanism supports smooth transitions between seated and standing positions.

Available in grey or beige, this lift chair is ideal for those seeking reliable comfort, local quality and functional ease.

Who it's for

- For clients at low risk of pressure injury who need a supportive, breathable seating solution.
- For therapists wanting a reliable, Tasmanian-made chair with motor and size options to suit a range of clients.

Key features

- **Available in single or twin motor** – allows for simultaneous or independent backrest and legrest control
- **Three size options** – supports a range of body types
- **Multi-layered pressure-reducing foam** – provides comfort and helps manage pressure care needs
- **Two-way stretch, breathable fabric** – anti-bacterial and waterproof to promote skin integrity
- **Easy-to-use hand control** – supports independence and comfort
- **Available in grey or beige** – neutral colours to suit different environments
- **Made in Tasmania** – high quality and locally produced

ADDITIONAL INFORMATION

Seat Height	490 mm
Brand	Topform Furniture
Seat Depth	470 mm
Seat Width	480 mm
Backrest Height	600 mm
Overall Width	780 mm
Overall Height	1020 mm
Max SWL	120 kg
User Height	153 cm to 160 cm
Motor	Single, Twin
Sizes	Mini, Medium, Maxi