



SKU:

TOPFORM LUXOR OPTIMA LIFTCHAIR





~~\$2,999.99~~ \$2,999.99

- Our Lift Chair Multi layer foam cushions developed specifically by health professionals to offer the optimum support and comfort for those requiring special pressure care.
- Quality Healthcare Pressure Management 2 way stretch fabric that offer numerous health benefits -anti bacterial, waterproof, breathable, anti microbial etc
- Topsuede denim on sides for a stylish and functional finish
- Available in Mini, Medium and Maxi size
- Available as 1 or 2 motor design.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Motor	Size
	S8330	\$2,999.99		Twin	Maxi
	S8327	\$2,799.99		Single	Maxi
	S8329	\$2,999.99		Twin	Medium
	S8328	\$2,999.99		Twin	Mini
	S8289	\$2,799.99		Single	Medium
	S8326	\$2,799.99		Single	Mini

PRODUCT DESCRIPTION

Ultrathane multi-way-stretch on contact surfaces which is:

- Extremely Hydrolysis Resistant
- Fluid Proof
- Anti-Bacterial
- Low Water Vapour Transmission
- Fire Retardant
- Hypoallergenic / Non-irritating
- Breathable
- Chemical Resistant for easy cleaning
- Available in Beige or Navy colours

Single vs twin vs quad motor lift chairs - what's the difference?

All ILS lift chairs have a lifting action, to support the user to stand from a seated position.

- Single motor chairs move the back rest and leg rest at the same time.
- Twin motor chairs allow the backrest and leg rest to move independently of each other, typically allowing for a greater degree of leg elevation and greater degree of recline compared to a standard single motor chair.
- Quad motor chairs allow independent movement of headrest, lumbar, backrest and footrest allowing for the greatest level of re-positioning of all chairs.

ADDITIONAL INFORMATION

Seat Height	490 mm
Seat Depth	470 mm
Seat Width	480 mm
Backrest Height	600 mm
Overall Width	780 mm
Overall Height	1020 mm
Max SWL	120 kg
User Height	153 cm to 160 cm
Weight capacity	100kg to 149kg
Motor	Twin, Single
Size	Maxi, Medium, Mini
Brand	Topform