




SKU: S13279

TRANSFER BENCH - SLIDING + SWIVEL BENCH ALUMINIUM RUST FREE

\$699.99

The Royale Medical Transfer Bench is a versatile transfer option, designed to improve quality of life and safety for users and attendees. The sturdy aluminium frame makes it a durable and lightweight showering solution with a safe working load of 130 kg. The transfer bench has height adjustable legs to alter seat height with durable anti-slip rubber tips. The seat, backrest and armrests are padded to provide ultimate comfort. The product comes with a sliding and swivel mechanism for the seat and has a drop-down feature for the armrests.



Scan QR code to
view product online

PRODUCT DESCRIPTION

The Royale Medical Transfer Bench allows users to remain seated whilst transferring in or out of the bathtub/shower. The product has a sturdy lightweight aluminium frame with a safe working load of 130 kg. The product has a comfortable padded blue seat, backrest, and armrests. The transfer bench has a swivel and sliding mechanism to assist users and attendees when transferring. The legs have adjustable height holes which will alter the seating height to meet user needs.

Key Features

- Lightweight aluminium frame
- Padded seat, backrest and armrests
- Swivel/Sliding seat
- Drop-down armrests
- Height adjustable seat
- SWL 130 kg
- Durable anti-slip rubber tips

ADDITIONAL INFORMATION

Brand	Royale Medical
Frame Material	Aluminium
Overall Width	84 cm
Overall Depth	52 cm
Seat Width	40 cm
Width Between the Arms	51 cm
Seat Height Range	50 cm – 60 cm
Armrest Height (from seat to top of arm)	12 cm
Weight Capacity/Max SWL	130 kg