




**SKU:** S13279

## **TRANSFER BENCH - SLIDING + SWIVEL BENCH ALUMINIUM RUST FREE**

**\$699.99**

The Royale Medical Transfer Bench is a versatile transfer option, designed to improve quality of life and safety for users and attendees. The sturdy aluminium frame makes it a durable and lightweight showering solution with a safe working load of 130 kg. The transfer bench has height adjustable legs to alter seat height with durable anti-slip rubber tips. The seat, backrest and armrests are padded to provide ultimate comfort. The product comes with a sliding and swivel mechanism for the seat and has a drop-down feature for the armrests.



Scan QR code to  
view product online

## **PRODUCT DESCRIPTION**

The Royale Medical Transfer Bench allows users to remain seated whilst transferring in or out of the bathtub/shower. The product has a sturdy lightweight aluminium frame with a safe working load of 130 kg. The product has a comfortable padded blue seat, backrest, and armrests. The transfer bench has a swivel and sliding mechanism to assist users and attendees when transferring. The legs have adjustable height holes which will alter the seating height to meet user needs.

### **Key Features**

- Lightweight aluminium frame
- Padded seat, backrest and armrests
- Swivel/Sliding seat
- Drop-down armrests
- Height adjustable seat
- SWL 130 kg
- Durable anti-slip rubber tips

## ADDITIONAL INFORMATION

<b>Brand</b>	Royale Medical
<b>Frame Material</b>	Aluminium
<b>Overall Width</b>	84 cm
<b>Overall Depth</b>	52 cm
<b>Seat Width</b>	40 cm
<b>Width Between the Arms</b>	51 cm
<b>Seat Height Range</b>	50 cm – 60 cm
<b>Armrest Height (from seat to top of arm)</b>	12 cm
<b>Weight Capacity/Max SWL</b>	130 kg