

SKU: S13279

ROYALE MEDICAL SLIDING & SWIVEL TRANSFER BENCH

\$699.99



The Royale Medical Sliding & Swivel Transfer Bench provides safe, controlled transfers for people with limited mobility who need assistance getting into and out of the bath. Its combination of sliding and swivel features minimises the need for lifting or twisting, enhancing safety and comfort.



Scan QR code to
view product online

GALLERY IMAGES





PRODUCT DESCRIPTION

The Royale Medical Sliding & Swivel Transfer Bench positions securely over the edge of the bathtub, enabling a person to sit and slide safely into the tub. The rotating seat allows easy positioning, while the sliding function supports smooth, controlled movement without the need to step over the bath edge – reducing the risk of slips, strains and falls.

The bench features a padded, contoured seat with backrest and armrests for comfort and stability. The seat locks securely at each end of the track for safer transfers.

Flip-up, lockable armrests offer added support during use and lift easily to enable unobstructed side transfers.

Height-adjustable legs accommodate individual needs and tub heights.

The rust-resistant frame ensures long-term durability, and suction-tipped feet enhance stability on wet bathroom floors.

Who it's for

- For clients who struggle with balance, strength or coordination and need a safer way to access the bath or shower.
- For therapists recommending a bench with swivel and slide functions to support safe, low-effort transfers.

Key features

- **Sliding and swivel seat** – simplifies bath entry with reduced effort
- **Padded seat, backrest and armrests** – enhance comfort and confidence
- **Flip-up, lockable armrests** – secure during use and lift for side transfers
 - **Lock-in-place mechanism** – stabilises the seat during transfer
- **Height-adjustable legs** – accommodate individual needs and tub heights
- **Non-slip suction feet** – provide stability on wet bathroom surfaces

- **Tool-free assembly** – quick and easy to set up or reposition

ADDITIONAL INFORMATION

Safe Working Load (SWL)	130 kg
Overall Width	84 cm
Overall Depth	52 cm
Seat Height	50 cm – 60 cm
Seat Width	40 cm
Width Between Arms	21 cm
Armrest Height	12 cm
Brand	Royale Medical

[Royale Medical Sliding & Swivel Transfer Bench User Manual](#)

Talk to us on the phone

Our experienced phone sales team can assist you with any product enquiries and help you make a purchase directly on the phone. We are available anytime between 9am and 5pm, Monday to Friday.

[1300 008 267](tel:1300008267)

No FAQ Found