



SKU:

ROYALE MEDICAL TRAVEL LITE POWER CHAIR

\$0.00

The Royale Travel Lite Folding Powered Wheelchair is a lightweight, foldable mobility solution for active lifestyles. Designed for easy travel and storage, it combines portability with a 150 kg weight capacity and up to 18 km range, promoting ongoing mobility and independence.



Scan QR code to
view product online

VARIATIONS

SKU	Price	Description	Sizes
S13647	\$3,999.99	The Royale Travel Lite Folding Powered Wheelchair is a lightweight, foldable mobility solution for active lifestyles. Designed for easy travel and storage, it combines portability with a 150 kg weight capacity and up to 18 km range, promoting ongoing mobility and independence.	8" Rear Wheels
S13648	\$4,499.99	The Royale Travel Lite Folding Powered Wheelchair is a lightweight, foldable mobility solution for active lifestyles. Designed for easy travel and storage, it combines portability with a 150 kg weight capacity and up to 18 km range, promoting ongoing mobility and independence.	10" Rear Wheels
S14225	\$4,999.99	The Royale Travel Lite Folding Powered Wheelchair is a lightweight, foldable mobility solution for active lifestyles. Designed for easy travel and storage, it combines portability with a 150 kg weight capacity and up to 18 km range, promoting ongoing mobility and independence.	12" Rear Wheels

PRODUCT DESCRIPTION

The Royale Travel Lite Folding Powered Wheelchair supports people who want to stay mobile and travel easily with a compact, motorised solution.

With a simple folding mechanism, the Travel Lite collapses in just a few steps for convenient storage in cars or tight spaces. Its lightweight frame and smart design make it well-suited to both everyday use and travel.

The chair offers a high safe working load of 150 kg, with a comfortable and adjustable seat cushion and backrest. Recline and armrest adjustments provide extra comfort.

Available in three rear castor sizes – 8", 10" and 12" – the Travel Lite is adaptable to different terrains and user preferences.

An intuitive joystick control enhances user confidence, and an optional attendant control is available for carers.

Solid, puncture-proof tyres reduce maintenance demands, while optional accessories such as width adjustment kits support customisation.

The Royale Travel Lite is airline-friendly, further supporting its use as a dependable travel companion.

Who it's for

- For clients seeking a portable powerchair for daily mobility and travel.
- For therapists wanting a folding powerchair with comfort, adjustability and transport ease for their client.

Key features

- **Lightweight folding design** – simplifies transport and storage
- **150 kg safe working load** – accommodates a broad range of users
- **Three castor size options** – to suit different needs and environments
 - **Adjustable backrest and armrests** – promote personal comfort
- **Ergonomic seat and backrest cushion** – for long-term seated support
- **User-friendly joystick control** – enhances confidence and independence
- **Optional attendant control** – supports shared use and carer involvement
 - **Solid, puncture-proof tyres** – for low-maintenance reliability
- **Width adjustment kits available as accessory** – to customise seating space
 - **Airline friendly** – suitable for travel

ADDITIONAL INFORMATION

Safe Working Load (SWL)	150 kg
Battery Type	Lithium
Battery Size	10 AH
Range	18 km
Max Speed	6.5 km/h
Size (W x L assembled)	65 cm x 100 cm
Seat Width	45 cm
Seat Depth	45 cm
Seat Height	52 cm
Weight	"33 kg (30.5 kg without batteries)"
Front Wheels	8" Castor Wheels
Rear Wheels	8", 10", 12" Options
Max Safe Incline Slope	9°
Turning Circle	120 cm
Brand	Royale Medical
Weight Of Each Rear Wheel	6.6 kg – 8" 7 kg – 10" 9 kg – 12"
Weight Of The Frame (without rear wheels)	17.4 kg (for all 3 sizes)
Weight Of Power Chair (without 24V10Ah battery)	30.6 kg – 8" 31.4 kg – 10" 35.4 kg – 12"
Total Weight Of Power Chair (with 24V10Ah battery)	33.1 kg – 8" 33.9 kg – 10" 37.9 kg – 12"
Unfolded Dimensions (LxWxH)	100 x 63 x 95 cm – 8" 100 x 65 x 98 cm – 10" 100 x 65 x 96 cm – 12"

Folded Dimensions (LxWxH)	80 × 63 × 35 cm – 8"
	80 × 65 × 35 cm – 10"
	80 × 65 × 35 cm – 12"
Seat Height (from ground)	52 cm
Arm Rest Height (from ground)	70 cm
Ground Clearence	10 cm
Wheel Base	47 cm
	45 cm
Wheel Tread	65 cm
Tyre Type	PU filled tyre
Back Rest Tilting Feature	5 inclinable angles
Folding/Unfolding Feature	1-second
Version	8", 10", 12"
Climbing Capability	9 Degrees

[Royale Medical Travel Lite User Manual](#)

[Royale Medical Travel Lite MSDS Certificate](#)

Talk to us on the phone

Our experienced phone sales team can assist you with any product enquiries and help you make a purchase directly on the phone. We are available anytime between 9am and 5pm, Monday to Friday.

[1300 008 267](tel:1300008267)

CAN THE WHEELS ON THE TRAVEL LITE ELECTRIC FOLDING POWER CHAIR BE A LARGER SIZE

Travel Lite Electric Folding Power Chair largest wheel size is 12.5 inches.