



SKU:

ROYALE MEDICAL TRAVEL LITE POWER CHAIR

\$0.00

The Royale Travel Lite Folding Powered Wheelchair is a lightweight, foldable mobility solution for active lifestyles. Designed for easy travel and storage, it combines portability with a 150 kg weight capacity and up to 18 km range, promoting ongoing mobility and independence.



Scan QR code to view product online



VARIATIONS

Image	SKU	Price	Description	Sizes
	S13647 \$	3,999.99		8" Rear Wheels
	S13648 \$	4,499.99		10" Rear Wheels
	S14225 \$	4,999.99		12" Rear Wheels



PRODUCT DESCRIPTION

The Royale Travel Lite Folding Powered Wheelchair supports people who want to stay mobile and travel easily with a compact, motorised solution.

With a simple folding mechanism, the Travel Lite collapses in just a few steps for convenient storage in cars or tight spaces. Its lightweight frame and smart design make it well-suited to both everyday use and travel.

The chair offers a high safe working load of 150 kg, with a comfortable and adjustable seat cushion and backrest. Recline and armrest adjustments provide extra comfort.

Available in three rear castor sizes – 8", 10" and 12" – the Travel Lite is adaptable to different terrains and user preferences.

An intuitive joystick control enhances user confidence, and an optional attendant control is available for carers.

Solid, puncture-proof tyres reduce maintenance demands, while optional accessories such as width adjustment kits support customisation.

The Royale Travel Lite is airline-friendly, further supporting its use as a dependable travel companion.

Who it's for

- For clients seeking a portable powerchair for daily mobility and travel.
- For therapists wanting a folding powerchair with comfort, adjustability and transport ease for their client.

Key features

- Lightweight folding design simplifies transport and storage
- 150 kg safe working load accommodates a broad range of users
- Three castor size options to suit different needs and environments
- Adjustable backrest and armrests promote personal comfort
- Ergonomic seat and backrest cushion for long-term seated support
- User-friendly joystick control enhances confidence and independence
- Optional attendant control supports shared use and carer involvement
- Solid, puncture-proof tyres for low-maintenance reliability
- Width adjustment kits available as accessory to customise seating space
- Airline friendly suitable for travel



ADDITIONAL INFORMATION

Safe Working Load (SWL) 150 kg

Battery Type Lithium

Battery Size 10 AH

Range 18 km

Max Speed 6.5 km/h

Size (W x L assembled) 65 cm x 100 cm

Seat Width 45 cm

Seat Depth 45 cm

Seat Height 52 cm

Weight "33 kg

(30.5 kg without batteries)"

Front Wheels 8" Castor Wheels

Rear Wheels 8", 10", 12" Options

Max Safe Incline Slope 9°

Turning Circle 120 cm

Brand Royale Medical

6.6 kg - 8"

Weight Of Each Rear Wheel 7 kg - 10"

9 kg - 12"

Weight Of The Frame

(without rear wheels)

17.4 kg (for all 3 sizes)



Waight Of Bawar Chair	30.6 kg - 8"
Weight Of Power Chair	31.4 kg - 10"
(without 24V10Ah battery)	35.4 kg - 12"
	33.1 kg - 8"

Total Weight Of Power Chair (with 24V10Ah battery) 33.1 kg - 8" 33.9 kg - 10" 37.9 kg - 12"

Unfolded Dimensions $100 \times 63 \times 95 \text{ cm} - 8"$ $100 \times 65 \times 98 \text{ cm} - 10"$ $100 \times 65 \times 96 \text{ cm} - 12"$

Seat Width 45 cm

Folded Dimensions (LxWxH) $80 \times 63 \times 35 \text{ cm} - 8"$ $80 \times 65 \times 35 \text{ cm} - 10"$ $80 \times 65 \times 35 \text{ cm} - 12"$

Seat Height (from ground) 52 cm

Arm Rest Height (from ground) 70 cm

Ground Clearence 10 cm
Wheel Base 47 cm
45 cm

Wheel Tread 65 cm

Tyre Type PU filled tyre

Back Rest Tilting Feature 5 inclinable angles

Folding/Unfolding Feature 1-second

Version 8", 10", 12"

Climbing Capability 9 Degrees