




SKU:

TRAVEL LITE ELECTRIC FOLDING POWER CHAIR - PHOSPHATE LITHIUM BATTERY




\$0.00

The Travel lite electric folding power wheelchair is a lightweight, easy folding mobility solution that will guarantee a comfortable ride. **The power wheelchair has a longer travel mileage and is designed to carry a load of 150 kg. Designed to fit inside any small space and a class-leading weight of only 17.4 kg.** Users can also enjoy a high-efficiency brushless hub motor and a joystick controller with 50% battery saving power. Its battery cycle is designed for 1000 charging cycles and is created out of pollution free materials.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Size
	S13647	\$3,999.99		8 inches
	S13648	\$4,499.99		10 Inches
	S14225	\$4,999.99		12 Inches

PRODUCT DESCRIPTION

Features

- Brushless Motors – Double Hub and Joystick Brushless. Latest market technology
- Huge travel range of up to 18 km with standard 10ah battery (included in unit) and 28 km with optional extra 15ah battery
- Huge 150 kg SWL design
- Tilt Recline – x5 multi adjust settings
- Smart Folding for easy storage
- Self-Adjustable Armrest Height and Angle
- Solid, Puncture-Proof Tyres
- Width adjustment kits available
- Airline Friendly and Safety Regulated. MSDS can be provided
- CE International Standards testing

ADDITIONAL INFORMATION

Brand	Royale Medical
Max Range	10 ah (18 km range) 15 ah (28 km range)
Max SWL	150 kg
Weight Of Each Rear Wheel	6.6 kg – 8” 7 kg – 10” 9 kg – 12”
Weight Of The Frame (without rear wheels)	17.4 kg (for all 3 sizes)
Weight Of Power Chair (without 24V10Ah battery)	30.6 kg – 8” 31.4 kg – 10” 35.4 kg – 12”
Total Weight Of Power Chair (with 24V10Ah battery)	33.1 kg – 8” 33.9 kg – 10” 37.9 kg – 12”
Unfolded Dimensions (LxWxH)	100 × 63 × 95 cm – 8” 100 × 65 × 98 cm – 10” 100 × 65 × 96 cm – 12”
Folded Dimensions (LxWxH)	80 × 63 × 35 cm – 8” 80 × 65 × 35 cm – 10” 80 × 65 × 35 cm – 12”
Seat Height (from ground)	52 cm
Arm Rest Height (from ground)	70 cm
Turning Radius	120 cm
Ground Clearance	10 cm

Front Wheels (Diameter*Width)	8" × 5 cm
Rear Wheels (Diameter*Width)	10" × 5 cm 12" × 5 cm
Wheel Base	47 cm 45 cm
Wheel Tread	65 cm
Tyre Type	PU filled tyre
Back Rest Tilting Feature	5 inclinable angles
Folding/Unfolding Feature	1-second
Max Speed	6.5 kmph
Version	8", 10", 12"
Climbing Capability	9 Degrees