




SKU:

ROYALE TRAVEL LITE FOLDING POWER CHAIR




\$0.00

The Royale Travel Lite folding powerchair is a lightweight, easy folding mobility solution that will guarantee a comfortable ride. Designed for people who want to stay active and mobile, the chair easily folds into a compact size for effortless storage and transporting purposes. The Travel Lite has a 18 km range (with standard battery) and SWL of 150 kg. Travel with ease, comfort, and independence with the Royale Travel Lite!



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Select Size
	S13647	\$3,999.99		8 inches
	S13648	\$4,499.99		10 Inches
	S14225	\$4,999.99		12 Inches

PRODUCT DESCRIPTION

The Royale Travel Lite is an ideal compact, lightweight, folding powerchair to assist user mobility and independence. The Travel Lite can be folded in just a few steps for convenient transport and storage purposes.

Built to last, the Royale Travel Lite has a high SWL of 150 kg with multiple standard adjustments that can be made and readily available accessories.

- Large foot platform to accommodate user needs
- 150 kg SWL
- Smart effortless folding design
- 3 different models available: 8", 10", 12" castors
- Ergonomic and comfortable seat cushion and backrest
- Backrest can be adjusted for the user to recline
- User friendly joystick
- Self-adjustable armrest height and angle
- Solid, puncture proof tyres
- Width adjustment kits available (as accessory)
- Attendee Control available (as accessory)
- Airline friendly

ADDITIONAL INFORMATION

Brand	Royale Medical
Max Range	10 ah (18 km range) 15 ah (28 km range)
Max SWL	150 kg
Weight Of Each Rear Wheel	6.6 kg – 8” 7 kg – 10” 9 kg – 12”
Weight Of The Frame (without rear wheels)	17.4 kg (for all 3 sizes)
Weight Of Power Chair (without 24V10Ah battery)	30.6 kg – 8” 31.4 kg – 10” 35.4 kg – 12”
Total Weight Of Power Chair (with 24V10Ah battery)	33.1 kg – 8” 33.9 kg – 10” 37.9 kg – 12”
Unfolded Dimensions (LxWxH)	100 × 63 × 95 cm – 8” 100 × 65 × 98 cm – 10” 100 × 65 × 96 cm – 12”
Folded Dimensions (LxWxH)	80 × 63 × 35 cm – 8” 80 × 65 × 35 cm – 10” 80 × 65 × 35 cm – 12”
Seat Height (from ground)	52 cm
Arm Rest Height (from ground)	70 cm
Turning Radius	120 cm
Ground Clearance	10 cm

Front Wheels (Diameter*Width)	8" × 5 cm
Rear Wheels (Diameter*Width)	10" × 5 cm 12" × 5 cm
Wheel Base	47 cm 45 cm
Wheel Tread	65 cm
Tyre Type	PU filled tyre
Back Rest Tilting Feature	5 inclinable angles
Folding/Unfolding Feature	1-second
Max Speed	6.5 kmph
Version	8", 10", 12"
Climbing Capability	9 Degrees