



TRAVEL LITE ELECTRIC FOLDING POWER CHAIR - PHOSPHATE LITHIUM BATTERY

\$2,999.99 - \$4,499.99

NEW & EXCLUSIVE TO ILS

The Travel lite electric folding power wheelchair is a lightweight, easy folding mobility solution that will guarantee a comfortable ride. **The power wheelchair has a longer travel mileage and is designed to carry a load of 150kg. Designed to fit inside any small space and a class-leading weight of only 15kg.** Users can also enjoy a high-efficiency brushless hub motor and a joystick controller with 50% battery saving power. Its battery cycle is designed for 1000 charging cycles and is created out of pollution free materials.

- Brushless Motors - Double Hub and Joystick Brushless. Latest market technology
- Phosphate Lithium batteries - platinum lifespan of more than 1000+ charges. 10 times the working life of standard lithium batteries.
 - Huge Travel range of up to 28km.
 - Huge 150kg SWL design
- Tilt Recline - x5 multi adjust settings.
 - Smart Folding for easy storage.
- Self-Adjustable Armrest Height and Angle
 - Solid, Puncture-Proof Tyres. Low Maintenance

- Compatible with our Hercules Lifter for Unit Transport if required.
 - Width adjustment kits available
- Airline Friendly and Safety Regulated.
MSDS can be provided
- CE International Standards testing



[How to Remove Wheels](#)




**FREE DELIVERY TO YOUR DOOR
AUSTRALIA WIDE (ONLINE ONLY)**



SKU: N/A

Category: [Power Chairs](#)

VARIATIONS

| Image | SKU | Price | Description | Size |
|---|--------|---|-------------|-------------|
|  | S13647 | \$4,495.99 <u>\$2,999.99</u> | | 8 inches |
|  | S13648 | \$4,995.99 <u>\$3,499.99</u> | | 10 inches |
|  | S14225 | \$5,495.99 <u>\$4,499.99</u> | | 12.5 inches |

GALLERY IMAGES



PRODUCT DESCRIPTION

ADDITIONAL INFORMATION

Max Range

10ah (18km range)

15ah (28km range)

| | |
|--------------------------------------|------------------------------------|
| Weight | 15kg (Rear wheels removed) |
| | 23kg - 27kg on total gross weight. |
| | 23kg - 8" |
| | 25kg - 10" |
| Max SWL | 27kg - 12" |
| | 150kg |
| | |
| Unfolded Dimensions (LxWxH) | 100*63*95cm - 8" |
| | 100*65*98cm - 10" |
| | 100*65*96cm - 12.5" |
| Folded Dimensions (LxWxH) | 80*63*35cm - 8" |
| | 80*65*35cm - 10" |
| | 80*65*35cm - 12.5" |
| Seat Height (from ground) | 52cm |
| Arm Rest Height (from ground) | 70cm |
| Turning Radius | 120cm |
| Ground Clearance | 10cm |
| Brand | Royale Medical |
| Front Wheels (Diameter*Width) | 8"*5cm |
| Rear Wheels (Diameter*Width) | 10"*5cm |
| | 12"*5cm |
| Wheel Base | 47cm |
| | 45cm |
| Wheel Tread | 65cm |
| Tyre Type | PU Filled tyre |
| Back Rest Tilting Feature | 5 inclinable angles |
| Folding/Unfolding Feature | 1-second |
| Max Speed | 6.5kmph |

Lifetime of brushless motor is 10 times of brush motor.
Lifetime of Phosphate Lithium batteries is 2~3 times of standard lithium batteries

Battery

Version

8", 10", 12.5"

Climbing Capability

9 Degrees