




SKU:

TROLLEY WALKER WITH HANDBRAKES - ADJUSTABLE HEIGHT

\$269.99

The Days adjustable trolley walker is supplied with a cream or white painted steel frame. Clip on trays can easily be removed for cleaning. Adjustable height easily accessible loop lockable brakes. As the trolley has brakes, this Trolley can also be used as a walking aid.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Colour
 A four-wheeled walker with a beige frame and two shelves. It has black handgrips and a front brake.	S1126	\$269.99		Beige
 A four-wheeled walker with a white frame and two shelves. It has black handgrips and a front brake.	S13000	\$269.99		White

PRODUCT DESCRIPTION

Finished in a cream or white painted steel frame, with handles that are adjustable in height the trolley walk with hand brakes is an essential daily accessory. It also has easily accessible loop lockable brakes.

Features

- Includes clip-on plastic trays that can be removed easily for cleaning.
- Can be used with or without trays or it can be used as a table.
- Rails on tray edges prevent items from slipping off.
- Lockable hand brakes mean this trolley can also be used as a walking aid and provides additional safety when standing from a seated position.
- Handles are adjustable in height.
- Easily accessible loop lockable brakes.
- Supplied with a cream painted steel frame.
- Comes in white or beige colour.

ADDITIONAL INFORMATION

Overall Width	58.5 cm
Overall Depth	56 cm
Handle Height Range	75.5 cm - 90.5 cm
Wheel Size	6"
Product Weight	8.25 kg
Maximum User Weight	150 kg