




**SKU:** S7827

## **TRUST CARE INDOOR WALKER WITH TRAY AND BAG**

**\$629.99**

The Trust Care Indoor Walker is a practical multi-functional rollator – you can use it as a walking aid, a kitchen trolley, or a table to eat from. It features clever features such as a removable plastic tray and one hand operated brakes. Easy to fold and transport, it is an ideal mobility aid for people who need extra support at home.



Scan QR code to  
view product online

## PRODUCT DESCRIPTION

The Trust Care Let's Go Indoor Walker is a lightweight easy-to-use rollator with exceptional function and safety features.

Key features include:

- Lightweight aluminium frame
- Adjustable height handles - The comfortable handles are covered with soft foam and are adjustable in height to suit your length.
- Manoeuvrability - Turns easily and can be used in narrow spaces like the bathroom and toilet.
- Multi-functional - In addition to being a walking aid, it can be used as a kitchen trolley, or a table to eat.
- Easy-to-use brake system - The brake is easy to operate with just one hand. The rollator brakes on both back wheels.
- Soft and quiet wheels - The wheels are made of soft rubber which make the rollator easy to turn, easy to roll, and use only a small space. The rollator moves easily over carpets.
- A removable plastic tray - Made from plastic with a matt surface, the tray can be placed tabletop height on the rollator. It is easy to lift off and on, and is used to transport a plate, cup of coffee or other smaller items.
- Textile basket - This may be placed under tray for storing small items.
- Easy to fold for transport or storage

## ADDITIONAL INFORMATION

|                      |                                        |
|----------------------|----------------------------------------|
| <b>Brand</b>         | Trust Care                             |
| <b>Handle height</b> | 78 cm – 90 cm                          |
| <b>Tray height</b>   | 68 cm                                  |
| <b>Length</b>        | 63 cm                                  |
| <b>Width</b>         | 52 cm                                  |
| <b>Weight</b>        | 6.5 kg                                 |
| <b>Max SWL</b>       | 100 kg                                 |
| <b>Frame type</b>    | Aluminium                              |
| <b>Colors</b>        | Beech/silver                           |
| <b>Warranty</b>      | 12 months, does not include wear parts |