




**SKU:** S1745

## **BATH RAIL - TUB MOUNT**

**\$149.99**

This powder coated steel bath rail is adjustable and clamps directly onto the side of most bathtubs. Suitable for users up to 100 kg and provides assistance for transferring in and out of the bath. A dual height grip position and clean aesthetic appearance round off the bath rail.



Scan QR code to  
view product online

## **PRODUCT DESCRIPTION**

A tri-wheel clamp ensures secure fastening and the rail is easy to install and suitable for use on most baths. Finished in a powder coated mild steel. Please note that this bath rail is only suitable for downward pressure.

## **ADDITIONAL INFORMATION**

<b>Adjustable Range</b>	75 mm - 175 mm
<b>Max SWL</b>	100 kg
<b>daily-living</b>	Bath Aids
<b>Brand</b>	Hero Medical