



SKU:

## HERO MEDICAL TUTOR PLUS FOREARM WALKER

~~\$999.99~~ \$999.99


The Hero Medical Tutor Plus Forearm Walker with 4 Wheels can very uniquely can be interchanged with rear legs skis. This means our customers – particularly within residential aged care and hospitals does not need to purchase another whole new unit fixed with a wheel/ski design for clients with different needs.

- This interchangeable option is exclusive to the HERO PLUS Forearm walker!
- A very economic option for streamlined equipment pools in hospitals and facilities.
- Designed to tutor the user into an upright gait.
- High 200kg Safe work Load – Extra Heavy Duty Sturdy Frame.
- Super Comfortable Durable Arm pads – Moulded Compound which allows Easy cleaning and Infection control. Great for facilities.
- Arm pads are studded into the framework ensuring pads remain in position at all times, unlike most makes on the market which consist of Velcro on Pads and Glue on Pads which do not remain in position for extended periods.



Scan QR code to  
view product online

**VARIATIONS**

Image	SKU	Price	Description	Size
 A blue metal walker with a black seat and handgrips, standard width.	S13823	\$549.99		Standard
 A blue metal walker with a black seat and handgrips, wider than the standard model.	S15395	\$999.99		Extra Wide

## PRODUCT DESCRIPTION

- Comes standard with x4 heavy duty 5" swivel castors - x2 standard at front and x1 rear locking castor / x1 rear directional castor. The directional castor allows a variety of different clinical benefits - particularly in a Rehab setting.
- Metallic Blue with a non-chip coated paint to allow for longevity in product condition
- Folds for easily for storage / transportation. Loops are featured on the folding hinges to allow a seamless folding process.
- Heavy duty Vinyl covering on the central frame to ensure the users fingers don't get caught and pinched during the folding motion of the frame.

### **Standard**

Total Height: 99.5-123.5cm

Total Width: 50.5cm

Gutter Width: 14cm

Gutter Length: 19cm

Handle Length: 37.5-47.5cm

SWL: 200kg

### **Extra Wide**

Total Height: 99.5-123.5cm

Total Width: 59.5cm

Gutter Width: 14cm

Gutter Length: 19cm

Handle Length: 37.5-47.5cm

SWL: 250kg

## **ADDITIONAL INFORMATION**

<b>Brand</b>	Hero Medical
<b>Max SWL</b>	200kg / Extra Wide SWL: 250kg
<b>Total height</b>	99.5-123.5cm / Extra Wide Total Height: 99.5-123.5cm
<b>Total width</b>	59.5cm / Extra Wide Total Width: 70cm
<b>Gutter width</b>	14cm / Extra Wide Gutter Width: 14cm
<b>Gutter length</b>	19cm / Extra Wide Gutter Length: 19cm
<b>Handle length</b>	37.5-47.5cm / Extra Wide Handle Length: 37.5-47.5cm
<b>Size</b>	Standard, Extra Wide