



SKU:

## HERO MEDICAL TUTOR PLUS FOREARM WALKER

~~\$999.99~~ **\$999.99**

The Hero Medical Tutor Plus Forearm Walker provides ergonomic upper body support for those unable to bear weight through their hands or wrists. Available in Standard and Extra Wide sizes, it supports safe ambulation with optimal posture.



Scan QR code to  
view product online

### VARIATIONS

| SKU    | Price    | Description | Sizes      |
|--------|----------|-------------|------------|
| S13823 | \$549.99 |             | Standard   |
| S15395 | \$999.99 |             | Extra Wide |

### PRODUCT DESCRIPTION

The Hero Medical Tutor Plus Forearm Walker is a height-adjustable walking aid designed for those who need forearm-based support. It helps reduce strain on the hands and wrists while encouraging upright, stable movement.

The walker includes moulded forearm platforms with padded supports and angle-adjustable handgrips for ergonomic control.

A powder-coated steel frame ensures durability, while swivel front castors, one rear locking castor and one rear directional castor enhances manoeuvrability.

Its design makes it suitable for both home and facility use, with easy-to-clean surfaces that support hygiene protocols and clinical safety.

### Who it's for

- For clients needing forearm-based weight-bearing support due to hand, wrist or lower limb limitations.
- For therapists prescribing an ergonomic walker that provides upper limb weight-bearing support.

### Key features

- **Padded forearm supports** – reduce wrist strain while allowing stable upper limb weight-bearing
  - **Height-adjustable** – customises fit for safety and postural alignment
  - **Angle-adjustable handgrips** – support ergonomic control during movement
- **Four-caster system** – includes two swivel front castors, one locking rear castor and one directional rear castor for manoeuvrability
  - **Powder-coated steel frame** – durable and suitable for clinical use
- **Supports upright walking posture** – improves balance and safety during ambulation

## ADDITIONAL INFORMATION

|                                |                      |
|--------------------------------|----------------------|
| <b>Sizes</b>                   | Standard, Extra Wide |
| <b>Safe Working Load (SWL)</b> | 200 kg               |
| <b>Overall Width</b>           | 50.5 cm, 59.5 cm     |
| <b>Overall Height</b>          | 99.5 cm – 123.5 cm   |
| <b>Gutter Width</b>            | 14 cm                |
| <b>Gutter Length</b>           | 19 cm                |
| <b>Handle length</b>           | 37.5 cm – 47.5 cm    |
| <b>Brand</b>                   | Hero Medical         |

[Hero Medical Tutor Plus Forearm Walker User Manual](#)

### **Talk to us on the phone**

Our experienced phone sales team can assist you with any product enquiries and help you make a purchase directly on the phone. We are available anytime between 9am and 5pm, Monday to Friday.

[1300 008 267](tel:1300008267)

No FAQ Found