



SKU:

HERO MEDICAL TUTOR PLUS FOREARM WALKER

~~\$999.99~~ \$999.99



The Hero Medical Tutor Plus Forearm Walker with 4 Wheels can very uniquely can be interchanged with rear legs skis. This means our customers – particularly within residential aged care and hospitals does not need to purchase another whole new unit fixed with a wheel/ski design for clients with different needs.

- This interchangeable option is exclusive to the HERO PLUS Forearm walker!
- A very economic option for streamlined equipment pools in hospitals and facilities.
- Designed to tutor the user into an upright gait.
- High 200kg Safe work Load – Extra Heavy Duty Sturdy Frame.
- Super Comfortable Durable Arm pads – Moulded Compound which allows Easy cleaning and Infection control. Great for facilities.
- Arm pads are studded into the framework ensuring pads remain in position at all times, unlike most makes on the market which consist of Velcro on Pads and Glue on Pads which do not remain in position for extended periods.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Size
 A blue metal walker with a black seat and handgrips, standard width.	S13823	\$549.99		Standard
 A blue metal walker with a black seat and handgrips, wider than the standard model.	S15395	\$999.99		Extra Wide

PRODUCT DESCRIPTION

- Comes standard with x4 heavy duty 5" swivel castors – x2 standard at front and x1 rear locking castor / x1 rear directional castor. The directional castor allows a variety of different clinical benefits – particularly in a Rehab setting.
- Metallic Blue with a non-chip coated paint to allow for longevity in product condition
- Folds for easily for storage / transportation. Loops are featured on the folding hinges to allow a seamless folding process.
- Heavy duty Vinyl covering on the central frame to ensure the users fingers don't get caught and pinched during the folding motion of the frame.

Standard

Total Height: 99.5-123.5cm

Total Width: 50.5cm

Gutter Width: 14cm

Gutter Length: 19cm

Handle Length: 37.5-47.5cm

SWL: 200kg

Extra Wide

Total Height: 99.5-123.5cm

Total Width: 59.5cm

Gutter Width: 14cm

Gutter Length: 19cm

Handle Length: 37.5-47.5cm

SWL: 250kg

ADDITIONAL INFORMATION

Brand	Hero Medical
Max SWL	200kg / Extra Wide SWL: 250kg
Total height	99.5-123.5cm / Extra Wide Total Height: 99.5-123.5cm
Total width	59.5cm / Extra Wide Total Width: 70cm
Gutter width	14cm / Extra Wide Gutter Width: 14cm
Gutter length	19cm / Extra Wide Gutter Length: 19cm
Handle length	37.5-47.5cm / Extra Wide Handle Length: 37.5-47.5cm
Size	Standard, Extra Wide