




**SKU:**

## **HERO MEDICAL ULTRA LIGHT ROLLATOR**

**\$199.99**

The Hero Ultra Light Rollator (walker) weighs 4.75kg and is designed using quality materials for greater comfort, safety and reliability. It has ergonomically designed hand grips, seat and backrest provides greater comfort to the user. The walker has multiple height adjustments and has a safe working load of 113kg to ensure that the walker is safe.



Scan QR code to  
view product online

VARIATIONS

| Image                                                                            | SKU    | Price    | Description | Colour |
|----------------------------------------------------------------------------------|--------|----------|-------------|--------|
|  | S14858 | \$199.99 |             | Gold   |
|  | S13781 | \$199.99 |             | Silver |

## **PRODUCT DESCRIPTION**

### **Hero Medical Ultra Light Rollator**

The Hero Ultra light Rollator is a streamlined walker that allows you to make the most of an independent life. The ergo form hand-grips are designed to give maximum comfort and safety without affecting the performance of the walker. The new comfort design of the padded seat and a backrest, the Hero Ultra Light Rollator provides comfort, safety and reliability all on one.

The Hero Ultralight Rollator's features include,

- Ergo-form handgrips fit more comfortably and place your hand at a natural angle for comfort and better brake access.
- Accessible brake system.
- Comfort Padded seat and backrest.
- Folds in position for easy storage or transport.
- Soft basket doesn't need to be removed when folding.
- Durable, ultra-light frame enables you to lift or move with minimal effort.
- 6" puncture free wheels are ideal for indoor and outdoor use.
- Reflective piping for added visibility and safety in low light situations.

## ADDITIONAL INFORMATION

|                                          |                                                                                       |
|------------------------------------------|---------------------------------------------------------------------------------------|
| <b>Brand</b>                             | Hero Medical                                                                          |
| <b>Weight Capacity / Max SWL</b>         | 113kg (249 lbs)                                                                       |
| <b>Weight without accessories</b>        | 4.75kg (10.4 lbs)                                                                     |
| <b>Overall height (floor to handles)</b> | Min 74cm – Max 92cm                                                                   |
| <b>Seating width</b>                     | 41cm                                                                                  |
| <b>Seat height (floor to seat)</b>       | Min 48cm – Max 58cm                                                                   |
| <b>Handle width</b>                      | 43cm(16.9 inches)                                                                     |
| <b>Length</b>                            | 79cm                                                                                  |
| <b>Castor size</b>                       | 6 inches                                                                              |
| <b>Colour</b>                            | Gold, Silver                                                                          |
| <b>Warranty</b>                          | 12 Months Warranty<br>Brake maintenance and adjustment is not covered under warranty. |