




**SKU:**

## HERO MEDICAL ULTRA LIGHT SEAT WALKER

**\$199.99**

The Hero Medical Ultra Light Seat Walker weighs just 4.75 kg and provides reliable, ergonomic support for everyday mobility. With padded seating, height-adjustable handles and a folding frame, it combines comfort, safety and convenience.



Scan QR code to  
view product online

### VARIATIONS

SKU	Price	Description	Colours
S14858	\$199.99		Gold
S13781	\$199.99		Silver

### GALLERY IMAGES



## PRODUCT DESCRIPTION

The Hero Medical Ultra Light Seat Walker is designed for individuals seeking dependable walking support without added weight. Built from high-quality materials, it delivers strong performance in a frame that's easy to handle and transport.

Ergo-form handgrips place the hands at a natural angle to reduce strain and improve access to the loop brake system.

The 6" puncture-free wheels support smooth indoor and outdoor use.

The padded seat and supportive backrest enhance comfort during rest breaks.

The walker folds upright for simple storage, with durable nylon bag in place.

### Who it's for

- For clients who need a lightweight, compact seat walker for stable, everyday mobility.
- For therapists prescribing an ultra-light walker that supports ergonomic posture and easy handling.

## Key features

- **Ergo-form handgrips** – support natural wrist positioning and easy brake access
  - **Accessible loop brakes** – ensure user control and responsive stopping
  - **Padded seat and backrest** – promote safe and comfortable rest
  - **Folding upright frame** – simplifies transport and storage
- **Nylon storage bag** – stays in place when folded for added convenience
  - **Durable ultra-light frame (4.75 kg)** – easy to lift and manoeuvre
- **6" puncture-free wheels** – support smooth indoor and outdoor movement

## ADDITIONAL INFORMATION

<b>Safe working load (SWL)</b>	113 kg
<b>Frame Weight</b>	4.75 kg
<b>Seat Height</b>	48 cm – 58 cm
<b>Seat Width</b>	41 cm
<b>Handle Height</b>	74 cm – 92 cm
<b>Front castor</b>	6"
<b>Rear castor</b>	6"
<b>Frame material</b>	Aluminium
<b>Colours</b>	Gold, Silver
<b>Brand</b>	Hero Medical

[Hero Medical Ultra Light Seat Walker User Manual](#)

### **Talk to us on the phone**

Our experienced phone sales team can assist you with any product enquiries and help you make a purchase directly on the phone. We are available anytime between 9am and 5pm, Monday to Friday.

[1300 008 267](tel:1300008267)