




SKU: B10036

DRIVE MEDICAL VIPER MOBILITY SCOOTER

From: \$4,999.99

The Viper Scooter offers modern styling and a whole host of excellent features. These include modern state of the art controls, one touch tiller adjustment and delta bars as standard. Safety has been the main concern during the development of this unique scooter.

The scooter features an anti-roll structure to avoid any difficulties accessing hills. A manual handbrake has been installed for emergency stopping power and an all-round powerful lighting package comes as standard including high-level indicators. **Twin 2A USB Ports for charging your accessories.**



Scan QR code to
view product online

PRODUCT DESCRIPTION

All New Viper Mobility Scooter - 75AH Platinum Gel Batteries

The Viper Scooter offers modern styling and a whole host of excellent features. These include modern state of the art controls, one touch tiller adjustment and delta bars as standard. Safety has been the main concern during the development of this unique scooter.

The scooter features an anti-roll structure to avoid any difficulties accessing hills. A manual handbrake has been installed for emergency stopping power and an all-round powerful lighting package comes as standard including high-level indicators. **Twin 2A USB Ports for charging your accessories.**

Key Features

- Top Speed of 12kph*
- Twin 2A USB Ports
- Maximum range of 48 kilometers
- All round suspension ensures a smooth ride and optimal handling
- 14" stylish alloy wheels with black pneumatic tyres for a smoother ride
- State of the art controls
- Comfortable rotating seat with sliding and reclining facility
- One touch hydraulic tiller adjustment to suit user requirements and assist easy transportation
- Delta bars as standard
- Height adjustable headrest to suit user requirements
- Armrests are angle adjustable and flip up for easy transfers
- Twin rear view mirrors for extra visibility
- Rear anti-tip wheels as standard
- Front removable storage basket as standard
- Freewheel facility allows movement of scooter without switching the motor on
- Large seat back storage pocket for storing essentials
- Additional hand operated cable brake for added safety
- Anti-roll back system safety device prevents scooter from rolling backwards on hills
- Batteries: 2 x 12v 75ah Platinum Gel Batteries
- Weight capacity 160kg

ADDITIONAL INFORMATION

Product Safe Working Load (SWL)	160kg (352.7 lbs)
Product Battery Size / Type	N/A (option of 75 AH Gel Batteries)
Product Range	45 km
Product Speed (Max)	10 km/hr
Product Colour	Red
Product Warranty	Mechanical 12 months Electrical 12 months Batteries 12 months
Product Size - in carton(s) (WxLxH) in cm	72 x 165 x 88
Product Size - in carton(s) (gross weight) in kg	170 kg
Product Size - overall assembled (WxLxH)* including pulled apart/folded for portable scooters in cm	56 x 143 x 132
Product Size - overall assembled (gross weight) in kg	120 kg
Product Size - seat base width in cm	47 cm
Product Size - seat base depth in cm	50 cm

Product Size - seat back width in cm	47 cm
Product Size - seat back height in cm	67 cm
Product Size - seat to floor in cm	62 cm – 64.5 cm
Product Size - front/ rear tyre type/size	14" / Pneumatic
Product Size - front/rear wheel size in inches	8"
Product Drive - turning circle in cm	134 cm
Product Drive - max safe ground clearance in cm	10 cm
Product Drive - max safe incline slope in degree	8°
Product Drive - suspension	Yes (Front and Rear)
Product Charger Type	5 AH
Product Motor Size	4 Pole / 800 watt
Product Lighting	LED
Product Display	Analog