



**SKU:** B10036

## **DRIVE MEDICAL VIPER MOBILITY SCOOTER**

**From: \$5,399.99**

The Viper scooter offers modern styling and a host of excellent features resulting in a smooth and comfortable driving experience. Just some of these great features include one touch infinite tiller adjustment and delta controls for easy steering and throttle control, a plush Captains Seat that is height adjustable and features a reclining backrest with adjustable headrest, 14" alloy wheels on front and rear suspension for a smooth ride on all terrain, extra bright LED headlights and turn signals, dual rear view mirrors, and a large front basket. For the ultimate driving experience don't settle for anything less than the Viper!



Scan QR code to  
view product online

## **PRODUCT DESCRIPTION**

- Front and rear suspension on 14" pneumatic tires for the smoothest ride possible
- One-touch infinite tiller adjustment
- Built-in USB port for charging your mobile devices
- Extra-bright LED lighting for superb visibility
- Ergonomic delta control handles for ease of steering and throttle control
- Large front plastic basket with carry handle
- Captain's seat is height adjustable and features a fold-down, reclining backrest with adjustable headrest
- Armrests are angle adjustable and flip-up for easy transfers
- LED headlights, parking lights, and turn signals for superior visibility
- Dual rear-view mirrors
- Anti-tip wheels
- Easy freewheel mode

## ADDITIONAL INFORMATION

<b>Product Safe Working Load (SWL)</b>	160kg (352.7 lbs)
<b>Product Battery Size / Type</b>	N/A ( option of 75 AH Gel Batteries)
<b>Product Range</b>	45 km
<b>Product Speed (Max)</b>	10 km/hr
<b>Product Colour</b>	Red
<b>Product Warranty</b>	Mechanical 12 months Electrical 12 months Batteries 12 months
<b>Product Size - in carton(s) (WxLxH) in cm</b>	72 x 165 x 88
<b>Product Size - in carton(s) (gross weight) in kg</b>	170 kg
<b>Product Size - overall assembled (WxLxH)* including pulled apart/folded for portable scooters in cm</b>	56 x 143 x 132
<b>Product Size - overall assembled (gross weight) in kg</b>	120 kg
<b>Product Size - seat base width in cm</b>	47 cm
<b>Product Size - seat base depth in cm</b>	50 cm

<b>Product Size - seat back width in cm</b>	47 cm
<b>Product Size - seat back height in cm</b>	67 cm
<b>Product Size - seat to floor in cm</b>	62 cm – 64.5 cm
<b>Product Size - front/ rear tyre type/size</b>	14" / Pneumatic
<b>Product Size - front/rear wheel size in inches</b>	8"
<b>Product Drive - turning circle in cm</b>	134 cm
<b>Product Drive - max safe ground clearance in cm</b>	10 cm
<b>Product Drive - max safe incline slope in degree</b>	8°
<b>Product Drive - suspension</b>	Yes (Front and Rear)
<b>Product Charger Type</b>	5 AH
<b>Product Motor Size</b>	4 Pole / 800 watt
<b>Product Lighting</b>	LED
<b>Product Display</b>	Analog