



SKU:

VULKAN/OPPO BOOT WALKER

\$99.99

Moon Boots are fracture boots designed to immobilise the foot and ankle following lower limb fractures, ligament strains and other injuries. They have two major advantages over standard casts.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Size
	S9763	\$99.99		Small
	S9764	\$99.99		Medium
	S9765	\$99.99		Large
	S9766	\$99.99		XLarge

PRODUCT DESCRIPTION

Since these boots use Velcro closures, they are can accomodate changes in swelling following an injury, as opposed to standard casts. Many people find lower limb casts become very loose in the final stages of use (due to disuse atrophy of muscles), but fracture boots can be progressively tightened to keep a snug fit.

Most fracture walkers fix the ankle at 90 degrees. However, some are articulated and the range of motion they allow at the ankle can be incrementally adjusted.

Patients using a cam walker are advised to wear a shoe of matching heel height on the opposite side.

The fracture boot heel is reasonably thick, and if you walk around lop-sided for several weeks it's common to develop some hip and back pain.

Product Features

- They can be taken off for showering
- They are designed to be walked in without the aid of crutches (however, depending on the injury, a patient may not be permitted to weight bear in the initial stages)

ADDITIONAL INFORMATION

Brand

Vulkan

Features:

- Lightweight and low profile aluminium side struts.
- Rocker-bottom design assists gait.
- Simple to apply and adjust by circumferential straps.