



SKU: S2168



PORTABLE WHEELCHAIR RAMP - EDGE BARRIER LIMITER

~~\$2,799.99~~ **\$2,799.99**

The Edge Barrier Limiter (E.B.L) range has the added security of edge barriers. Available in different lengths to provide solutions to obstacles from 22.5 cm to 87.5 cm in height. Particularly suited to providing instant access to buses, trains and buildings where a permanent ramp is impractical. A simple concertina folding action enables the Invacare Portable Ramp to be deployed in seconds.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Length in cm
	S2168	\$799.99		090
	S2209	\$999.99		120
	S2193	\$1,049.99		135
	S7648	\$1,149.99		150
	S7547	\$1,199.99		165
	S0329	\$1,699.99		200
	S2262	\$2,299.99		250
	S2255	\$2,799.99		300

PRODUCT DESCRIPTION

- **Portable** – Invacare Portable Ramps can be quickly folded into a single compact unit – allowing for easy storage or transport in any vehicle or home environment.
- **Lightweight** – Made from a unique combination of unidirectional glass and carbon fibre materials. Invacare Portable Ramps are very strong, yet light enough to be portable.
- **Compact** – Multi-paneled design held together with a patented hinge system connected by webbing enabling Invacare Portable Ramps to fold with a concertina action, allowing for a much smaller width when folded.
- **Safe** – Clearly marked and grip tape covered, Invacare Portable Ramps are also fitted with pads at the top and base. Invacare Portable Ramps provide a non-slip surface, for all weather conditions.
- **Features**- 5 cm Edge barrier on each side.

Ramp Sizes

Length	Width	Weight
90cm	82cm	6kg
120cm	82cm	9.7kg
135cm	82cm	NA
150cm	82cm	NA
165cm	82cm	13kg
200cm	82cm	13.7kg
250cm	82cm	NA
300cm	82cm	27.3kg

ADDITIONAL INFORMATION

Brand	Invacare
Length in cm	090, 120, 135, 150, 165, 200, 250, 300
Ramp width	82 cm
Safe working load	300kg