




SKU:

RAMP - INVACARE PORTABLE DECPAC EBL

~~\$3,699.99~~ **\$3,699.99**

A ramp for those who like to get about choose this EBL model. The Decpac EBL range has been developed to give users access to shops, offices and the community. It has a 5 cm edge barrier design for safety and security to prevent steering off the edge of the ramp. These ramps provide instant access into buses, trains or other facilities via the concertina folding design that is deployed in seconds!

Invacare Portable Ramps are manufactured with a unique combination of uni-directional glass and carbon fibre materials to create an extremely strong and lightweight ramp suitable for a wide range of uses.



Scan QR code to
view product online

VARIATIONS





Image	SKU	Price	Description	Select size
	S2168	\$849.99		90 cm
	S2209	\$1,049.99		120 cm
	S2193	\$1,099.99		135 cm
	S7648	\$1,149.99		150 cm





Image	SKU	Price	Description	Select size
	S7547	\$1,249.99		165 cm
	S0329	\$1,899.99		200 cm
	S2262	\$2,499.99		250 cm
	S2255	\$2,999.99		300 cm

Image	SKU	Price	Description	Select size
	S8068	\$3,699.99		350 cm

PRODUCT DESCRIPTION

- EBL - Edge Barrier Limiter feature is fitted to these ramps to provide additional safety
- Portable - Invacare Portable Ramps can be quickly folded into a single compact unit, allowing for easy storage or transport in any vehicle or home environment
- Lightweight - Made from a unique combination of unidirectional glass and carbon fibre materials. Invacare Portable Ramps are very strong, yet light enough to be portable
- Compact - Multi-panelled design held together with a patented hinge system connected by webbing enabling Invacare Portable Ramps to fold with a concertina action, allowing for a much smaller width when folded
- Safe - Clearly marked and grip tape covered, Invacare Portable Ramps are also fitted with pads at the top and base. Invacare Portable Ramps provide a non-slip surface for all weather conditions
- Strong - Invacare Portable Ramps are tested to AS/NZ Standards, up to 300 kg for security and peace of mind
- Simple to use - Invacare Portable Ramps require no assembly allowing the ramp to be deployed in seconds, and immediately ready for use!

ADDITIONAL INFORMATION

Brand	Invacare
Ramp width	82 cm
Safe working load	300kg