




SKU:

WONDER SHEET PLUS FITTED SHEET

~~\$160.00~~ \$160.00

The Wonder Sheet Plus Fitted Sheet combines a satin centre panel with cotton/polyester sides for enhanced bed mobility. Its design leaves the mattress edges fabric-only, providing added stability when moving near the sides of the bed.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Depth
	S314004	\$119.99	Single 35 cm	
	S314415	\$119.99	Long Single 35 cm	
	S313707	\$129.99	King Single 35 cm	
	S324253	\$149.99	Double 35 cm	
	S323457	\$149.99	Long Double 35 cm	
	S313147	\$160.00	Queen 45 cm	

PRODUCT DESCRIPTION

The Wonder Sheet Plus Fitted Sheet is a therapeutic bedding solution that helps people move more easily in bed while offering added support near the edges. Its 90 cm satin panel is centred under the shoulders, back and hips, but stops 17 cm short of each side of the bed. This makes it easier to reposition and pivot, while retaining grip and stability along the mattress edges.

This design is particularly helpful for people who need smooth sliding in the centre, but extra traction when getting in and out of bed. It reduces the effort and discomfort involved in turning or adjusting position, especially for those with reduced mobility or pain.

Made from a glossy charmeuse satin and 50/50 cotton/polyester blend, it's both soft and durable. The sheet is available in two fitted depths to suit a range of mattress sizes.

Who it's for

- For clients who want friction-free movement in bed with added edge stability.
- For therapists seeking a repositioning aid that still supports safe bed entry and exit.

Key features

- **Satin panel with fabric edges** – supports in-bed mobility while maintaining edge stability
- **90 cm long satin section** – placed under the torso and hips to ease movement
- **Glossy charmeuse satin and cotton/polyester blend** – combines comfort and durability
- **Fitted sheet with depth options** – 35 cm depth for mattresses up to 30 cm; 45 cm depth for mattresses 30-40 cm

ADDITIONAL INFORMATION

Brand

Neeki Designs