


SKU:

XENON2 DUAL

Whether you're looking for one of the very best in driving dynamics or compact folding dimensions, there's a Xenon2 just for you.

Will it be the clean, streamlined design of the Xenon2 FF? Alternatively how about the reinforced rigidity of the Xenon2 Dual? Or the practical portability of the Xenon2 SA? Whichever one of the three feature-packed models you choose, the sophisticated elegance of the Xenon2 will always make you the centre of attention.

Ultra lightweight wheelchairs with outstanding folding and driving performance.



Scan QR code to
view product online



PRODUCT DESCRIPTION

| | |
|---------------------------------|-----------------------------------|
| Max. user weight: | 140 kg |
| Seat width: | 32 – 50 cm |
| Seat depth: | 34 – 50 cm |
| Seat height front: | 43 – 55 cm |
| Seat height rear: | 37 – 50 cm |
| Backrest height: | 25 – 47.5 cm |
| Backrest angle: | 75° to 103° in 4° increments |
| Centre of gravity (COG): | 4 – 13 cm |
| Camber: | 0° / 2° / 4° |
| Frame Version: | Reinforced frame with fixed front |
| Chair weight: | from 9.5 kg |
| Frame colours: | 32 |
| Crash tested: | yes |

ADDITIONAL INFORMATION

Brand

Sunrise Medical