

# Jazzy Air<sup>2</sup>



Shown in Pearl White

**136 kg**

Weight Capacity

Up To

**6.4 km/h**

Maximum Speed

Elevated & Non-Elevated

**61.5 cm**

Turning Radius

Up To

**27 km**

Range Per Charge\*



Bright LED marker lights

At the touch of a switch, the **Jazzy Air<sup>2</sup>** elevates 30.4 cm in just 11 seconds. With a sleek, more modern look, this power chair safely drives at speeds up to 6.4 km/h while elevated, allowing users to enjoy face-to-face social engagement. With Mid-Wheel 6<sup>®</sup> Drive and Active-Trac<sup>®</sup> Suspension, the Jazzy Air 2 offers superb performance and stability both indoors and out.

## Specifications

Model Name	Jazzy Air® 2
Weight Capacity	136 kg (300 lbs.)
Maximum Speed <sup>1</sup>	6.4 km/h (4 mph)
Ground Clearance <sup>2</sup>	2.54 cm (1") at motor bracket
Turning Radius <sup>1</sup>	61.5 cm (24.25")
Overall Length <sup>2</sup>	106.7 cm (42")
Overall Width <sup>2</sup>	64.1 cm (25.25")
Front Wheels	15.24 cm (6") flat free
Drive Wheels	25.4 cm (10") flat free
Rear Wheels	15.24 cm (6") flat free
Drivetrain	Two motor, Mid-Wheel 6®
Braking System	Regenerative and electromechanical
Suspension Type	Active-Trac (ATX) with Mid-Wheel 6® Technology
High-Back Seat Size	45.7 x 45.7 - 50.8 cm (18" x 18" - 20")
Seat-to-Floor Range	54 cm - 83.8 cm (21.25" - 33")
Range Per Charge <sup>2,4</sup>	Up to 27 km (16.8 miles)
Standard Electronics	40A, PG VR2
Component Weights <sup>6</sup>	Base: 40 AH Batteries 120 kg 40 AH Batteries: 14.28 kg
Battery Requirements <sup>3,5</sup>	(2) 12V, 40 AH
Battery Charger	Off-board, 3.5 amp
Warranty	Frame: Lifetime limited Drivetrain: 2-year Electronics: 1-year Battery: (Gel battery) 1-year (AGM battery) 6-months

1. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of (+ or -) 10%.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
3. AGM or gel cell type required.
4. Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
5. Battery weight may vary based on battery manufacturer.
6. Weight includes base, seat, and batteries.

The information contained herein is correct at time of publication; we reserve the right to alter specifications without prior notice.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition, and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius, and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

## Elevate your life to the next level!

- Newly re-designed lift features 30.4 cm of power adjustable seat height in 11 seconds
- Convenient "Air" switch on armrest activates elevation
- Active-Trac® Suspension with Mid-Wheel 6® Technology allows for maximum stability in the elevated and non-elevated positions
- Bright LED marker lights offers visibility in all environments
- High-back or solid seat pan, reclining sport seat
- Height and angle adjustable headrest
- Front under-seat storage keeps your valuables accessible

## Additional Features

- Angle, width, height and depth adjustable armrests
- Integrated foot platform with rubber insert and 5 cm of adjustment for comfort in both the elevated and non-elevated positions
- Single-sided castor forks for enhanced manoeuvrability

## Available Accessories

- Rear bag / walking stick holder
- Walker holder\*
- Oxygen tank holder\*
- Mobile phone holder
- Cup holder
- Saddle Bag
- Under-seat storage compartment

\*Must have Pride rear bag to attach

## Available Colours

