

Q6 EDGE 3 *Stretto*™



Q6 EDGE 3 STRETTO™
shown with iLevel® power seat elevation,
Q-Logic 3 environmental controls, and
Stealth® TRU-Comfort 2 cushion.



The Q6 Edge 3 Stretto™ features a narrow width for outstanding maneuverability and is equipped with independent SRS (Smooth Ride Suspension), making it ideal for children, teenagers and small adults. The Stretto is available with 31.8 cm or 35.6 cm drive wheels and optional 30.5 cm of iLevel® power adjustable seat height.

Narrowest
complex
rehab base
on the market




QUANTUM®
#1 FOR REHAB POWER



quantumrehab.com

Specifications

DRIVE WHEELS	35.6 cm or 31.8 cm pneumatic or flat free tires
CASTER WHEELS	Front: 12.7 cm flat free Rear: 15.2 cm flat free
MAXIMUM SPEED ¹	Up to 10 km/h Up to 5.6 km/h w/iLevel activated
GROUND CLEARANCE	7.4 cm
TURNING RADIUS ²	53.3 cm
OVERALL LENGTH ²	82.5 cm without front riggings
BASE WIDTH	55.25 cm with 35.6 cm drive wheels 52.71 cm with 31.8 cm drive wheels
SEATING SIZES	
TRU-Balance® 3 Static Seat	Width: 30.5 cm-51 cm Depth: 30.5 cm-51 cm
TRU-Balance® 3 Power Tilt	Width: 30.5 cm-51 cm Depth: 30.5 cm-51 cm
TRU-Balance® 3 Tilt and iLevel®	Width: 30.5 cm-51 cm Depth: 30.5 cm-51 cm
SEAT-TO-FLOOR HEIGHTS	
TRU-Balance® 3 Static Seat ⁵	40.6 cm - 45.7 cm
TRU-Balance® 3 Power Tilt ⁵	44.5 cm - 49.5 cm
TRU-Balance® 3 Tilt & iLevel® ⁵	44.5 cm - 49.5 cm
MANUAL RECLINE ⁶	Yes (TRU-Balance® Flex only)
BATTERY SIZE ³	40 AH (2) or 55 AH (2)
BATTERY WEIGHT ⁷	40 AH: 14 kg each 55 AH: 18.1 kg each
AVAILABLE ELECTRONICS ⁸	120A Q-Logic 3 EX 75A Q-Logic NE 75A Q-Logic NE+
BATTERY CHARGER	8A Off-board
MOTOR PACKAGES	2-Pole, Mid-Wheel ⁶
WEIGHT CAPACITY ⁴	136 kg
BASE WEIGHT WITHOUT SEATING SYSTEM	57.3 kg

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
- 5) All measurements +/- 0.3 cm with 45.7 cm x 45.7 cm seat, fixed towers, 136 kg weight capacity and solid tires (add 0.95 cm with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.
- 7) Battery weight may vary +/- 1.4 kg based on manufacturer.
- 8) All Quantum Rehab® units can be ordered with specialty controls.
- 9) Call for availability.

Features

- Independent SRS (Smooth Ride Suspension)
- iLevel® option available with up to 30.5 cm of seat elevation at walking speed (up to 5.6 km/h)
- Fender lights and USB charger included
- Made in America of foreign and domestic parts

Options

- iLevel® option allows 30.5 cm of seat elevation while at walking speed (5.6 km/h)
- Q-Logic 3 EX Drive Controls
- TRU-Balance® 3 Seating
- Trilogy vent tray⁹
- Compact vent tray⁹
- 30.5 cm x 30.5 cm seat sizes with iLevel (SOF only)

Colours*

