



BC2 Super Kingsize

Specifications

Overall Width: 94 cm

Overall Depth: 74 cm

Overall Height: 99—114 cm

Seat Width: 74 cm

Seat Depth: 47 cm

Seat Height: 41 —56 cm

Back Height: 72 cm

SWL: 300 kgs

Warranty: 10 years structural

A bariatric chair featuring a fully welded steel frame with adjustable leg height, fire retardant foams, 3 segmented back cushion and a one piece ergonomically designed shell. Back and seat cushion are loose for ease of cleaning.

Optional extras

- rear transport wheels
- Side wings