

The Bellavita bathlift reclines further, travels lower and is the lightest in the Drive DeVilbiss range of battery-powered reclining bathlifts\* It is easy to install, simple to operate and offers comfort, safety, reliability and convenience.





## **Features**

- At the touch of a button, gently lowers users down to bottom of bath for bathing and relaxation then raises them to top again for safely getting out
- The backrest reclines to a maximum angle of 50 degrees, operated easily by the hand control and enabling the user to relax in comfort
- Front aperture in seat improves access for personal hygiene
- Can be used with or without covers; washable hygienic covers are easily fitted/removed and provide comfort when bathing (available in white or blue)
- Very low seat travel of 6cm (2.4") from bottom of bath allows improved water immersion, ideal for shallow bath
- A seat height of 48cm (19") enabling easy and safe transfers in and out of the bath

- At 9.3kg\*\* the Bellavita is very lightweight making it easy to install and remove from the bath
- Separates easily and quickly into 2 lightweight sections for ease of fitting/removal, storage and transportation
- Contoured, slim-line profile further facilitates handling and gives a modern, aesthetically pleasing appearance to suit most bathrooms
- 4 quick release suction feet secure the lift in the bath. The front 2 are self-venting and flexible enabling easier install and added suction to the bath
- Lightweight, fully waterproof and floatable hand control with light touch buttons is easily operated

- Hand control recharges quickly for convenience and multiple use environments
- Indicator lights signal when charging is required and built in safety feature prevents bathlift from lowering without sufficient charge to return it to top of bath
- Emergency stop button is incorporated into hand control for safety
- Lightweight plastic frame is easy to clean
- Easily fits into most baths
- Wide range of optional extras, including a carry bag, covers, side flap sizes, belt systems and side flap bumpers

## **Specifications**

Overall Length: 83cm (33")

Overall Width: 70cm (28")

Seat Width: 35cm (14")

Seat Depth: 42cm (17")

Overall Seat Length: 57cm (22")

Sitting Height Min: 6cm (2")

Sitting Height Max: 48cm (19")

Backrest Height: 64cm (25")

Backrest Width: 32-36cm (13-14")

Backrest Reclines: 50°

Baseplate Width: 31cm (12")

Baseplate Length: 56cm (22")

Load Capacity Max: 140kg (22st)

Total Weight

(excluding handset): 9.3kg (20.5lb)

Lightweight, fully waterproof floating hand control

Backrest reclines to an incredible 50°

**Leaflet Code: LL127** 

If you would like to view this product please contact your local stockist.

We have a continuous product improvement policy and consequently reserve the right to amend design and specification without prior notice. All sizes and weights stated are nominal. All errors and omissions excepted. Due to the limitations of photography and the printing process colours shown may not be 100% accurate. \*\*excluding handset. REV4/13022020

