

 Breezy



BREEZY EC WHEELCHAIRS

The best combination of durability, selection and value.



Breezy EC 2000

A quality wheelchair with options offering the flexibility to accommodate a variety of users.

Standard Features

- Dual axle allows the seat to be easily lowered to hemi height
- Chrome-plated, carbon steel frame is built to last, great for rental fleets
- Composite mag wheels are molded with contours to enhance strength and durability
- Full size chart pocket
- Industry-compatible front riggings
- 113.5 kgs weight capacity (136 kgs on 20" seat width)

Options for Versatility

- Removable desk-length
- Swing-away footrests

Models

- 216RADPS 16" width, removable padded desk arms, swing-away footrests
218RADPS 18" width, removable padded desk arms, swing-away footrests
220RADPS 20" width, removable padded desk arms, swing-away footrests



Breezy EC 2000HD Heavy-Duty

Extra wide seat widths and reinforced construction accommodate special needs.

Features

- Heavy-duty reinforced chrome steel frame supports weights up to 204.5 kgs
- Seat widths of 22" or 24" accommodate larger users.
- Durable, double-embossed vinyl upholstery
- Swing-away footrests
- Oversized chart pocket
- Removable desk length, armrests

Models

- 2HD22RADPS 22" width, padded desk arms, swing-away footrests
2HD24RADPS 24" width, padded desk arms, swing-away footrests





Breezy EC 2000HD Bariatric

Extra roomy seat widths and rugged steel construction provide mobility for users up to 318 kgs

Features

- Available in 26", or 28" seat width with 20" seat depth. supports weights up to 318 kgs.
- Rugged construction with reinforced cross frame supports weight: up to 318 kgs
- Detachable padded desk arm can be reversed allowing easier access for propelling
- Industry-compatible front riggings
- Durable nylon upholstery
- Aluminum swing-away footrests

Models

2HD26RADPS	26" width, padded desk arms, swing-away footrests
2HD28RADPS	28" width, padded desk arms, swing-away footrests



JAY® Seating Options

JAY® Basic Back

The Jay Basic Back has been redesigned to easily adapt to standard wheelchairs. Choose strap-style hardware for easy attachment - just strap around wheelchair back canes. The clamp-style hardware provides 0° and 8° of angle adjustability and snaps easily into place.


- Designed for clients with mild postural needs
- Provides mild posterior pelvic-stabilization and thoracic support
- Three height options provide a more personalized fit
- Lightweight for improved propelling efficiency and easy transport
- Lightweight, contoured, aluminum shell and soft foam insert



JAY® Basic Cushion

The Jay Basic provides a soft and durable support for clients with mild postural needs who are at low risk for skin breakdown.

- Beveled to fit between the seat rails on wheelchairs with seat upholstery
- Durable molded foam provides mild contouring for increased sitting tolerance and postural stability
- Provides stability to reduce sliding and shearing
- Available for chair widths 14" to 24"



SPECIFICATIONS

	EC 2000	EC 2000HD	EC 2000HD Bariatric
Seat width	16"/18"/20"	22"/24"	26"/28"
Seat Depth	16"	18"/18"	20"
Overall Width			
Fixed Arm	22.5"/24.5"	N/A	N/A
Removable Arm	24.5"/26.5"/29.5"	30.5"/32.5"	35"/37"/39"
Overall Depth			
Fixed Arm, Fixed F/R	41.5"	N/A	N/A
Fixed Arm (w/o riggings)	32"	N/A	N/A
Removable (w/o riggings)	29"	29"/32"/32"	37"
Seat to Floor Height	19.75"	19.75"	19.75"
Back Style	Fixed	Fixed	Fixed
Back Height	16"	16"	18"
Arm Style	Removable Desk	Removable Desk	Removable Desk
Front Riggings	Fixed Swing-away	Swing-away	Swing-away Elevating
Rear Axle	Dual Position	Single Position	Single Position
Rear Wheels	Composite Mag	Composite Mag	Steel Mag
Hand Rims	Composite	Chromed Steel	Chromed Steel
Upholstery	Vinyl	Vinyl	Nylon
Product Weight (w/o riggings)	18.5 kgs	26 kgs	36 kgs
Weight Capacity	113 / 136 kgs	204.5 kgs	318 kgs
Warranty			
Frame	3 Years	3 Years	3 Years
Tires	None	None	None