

Softform Premier

DURABILITY AND PROTECTION



Equipped with a suite of features that ensures durability, sustained high degree of protection against tissue damage for the patient and an ingeniously designed insert pad, the Softform Premier mattress really does provide exceptional pressure reduction and comfort.

High-Very High Risk



Yes, you can.®

Making Life's Experiences Possible



Softform Premier

Unique Features

Benefits

Deep contoured insert pad for enhanced pressure care area	Maximises body contact area and optimises pressure reduction
Cover fabric conforms to deep contoured foam castellations on the upper surface of the insert pad when mattress is profiled	Minimizes rucking of cover fabric on top surface, reducing potential for tissue damage
Suitable for High Risk patients / Very High Risk when used in conjunction with a profiling bed	Can contribute significantly to pressure ulcer prevention
No need to turn over	Reduces manual handling issues
Single piece foam 'U' core	Facilitates patient transfer by providing effective side-wall support
Double-coated poly-urethane sides	Reduces the risk of strike through, hydrolysis or delamination
Two way stretch polyurethane waterproof cover	Provides excellent recovery and durability and helps reduce shear and friction forces
High frequency welded seams and zip	Reduces risk of fluid ingress
Full length waterfall flap conceals the zip	Further risk reduction of fluid ingress (through zip)
Toughened PU coated base	Prolongs longevity of the mattress
Water penetration resistance	35 kPa minimum (BS 3424: part 26). 35 kPa minimum after 20,000 cycles Martindale Abrasion Test EN ISO 5470-2

Tech Specs

Codes

Single	User weight 247.6kg	Length 1970mm	Width 880mm	Depth 152mm
King Single	247.6kg	1970mm	1070mm	152mm

WARRANTY

Warranty
4yr cover
8yr foam

- Softform Premier Single - SPM197/88
- Softform Premier King Single - SPM197/107



FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz

Australia: Freephone 1800 460 460, email orders@invacare.com.au or visit www.invacare.com.au

Asia: Phone +61 2 8839 5330, email anzexportsales@invacare.com or visit www.invacare.com

Thailand: Phone +66 0 2 821 5515 or visit www.invacare.co.th



Yes, you can.®