


Icon Plus



Icon Plus

Pressure Care Support Surface



Acute Care Medium Risk Static Care Support Surface

The Icon Plus provides true leadership in static pressure care in an acute preventative setting. Tried and trusted within leading hospital facilities nation-wide, the Icon Plus offers outstanding durability, and an unrivaled consistency of performance. The Icon Plus is hand made in Australia and offers superior patient comfort, effective pressure relief and skin protection for patient up to medium risk. An extensive variety of configurations are available providing for a wide range of patient needs and requirements including custom sizing for precision fitting to the bed base used within the care setting.



APPLICATIONS

- ✓ Hospital Ward Usage
- ✓ Aged Care Facilities
- ✓ Home Care




AIRO-FORM IMMERSION UPPER LAYER

- ✓ Alleviates localised pressure points
- ✓ Provides maximum contact with patient for distribution of interface pressure
- ✓ High temperature stability to deliver consistent immersion, performance and comfort
- ✓ Unique Airo-form ventilation system minimizes heat entrapment and enhances core circulation providing a cooler sleep surface
- ✓ Visco-elastic cell structure provides low-shear surface




PROFILE STABILISING SIDES

- ✓ Enables ease of primary patient care
- ✓ Facilitates patient transfers
- ✓ Unique hinge points allow mattress to easily conform to varied bed positions
- ✓ Encourages central patient positioning



PREMIUM GRADE MID-LAYER

- ✓ Open cell high quality foam
- ✓ Provides exceptional comfort and postural support for repositioning and patient care



PREMIUM FIRM BASE LAYER

- ✓ Provides protection against 'bottoming-out' during transfer, repositioning and whilst in the inclined/fowler position
- ✓ Provides exceptional comfort, support for repositioning and patient care

TRIPLE LAYERED CORE


- ✓ Foam conforms to AS/NZS 4088.1 (Int):1993 and FMVSS302 & SAE T2002
- ✓ Anti-microbial treatment applied to foam maximises longevity and core micro-climate
- ✓ Conforms to EPUAP, NPUAP & PPIA Guidelines for quality and comfort

Icon Plus

SIZE	CODE
1980x880x150	S16398
2030x880x150	S16399
2030x1050x150	S16400
1980x1050x150	S16401


Effective Weight Capacity: 40-165KG
Cover Type: Advantiflex Classic Cover

Pressure Relief




7/10

Shear Protection



5/10

Breathability



6/10

QUICK GUIDE


Forté Pressure Care Covers

Superior Quality, High Performance Medical Grade Pressure Relieving Cover



Advantiflex Classic Icon Plus Cover

Belgium made polyurethane coated upper and heavy duty PVC Base ensures that Forté Advantiflex Classic Covers deliver consistent care, comfort and durability.



Premiflex Classic Sovereign Plus Cover

Premiflex Classic Covers are differentiated by their superior performance materials. This includes heavier weight fabrics, increased shear protection, greater resistance to cleaning agents and superior durability within the acute care setting.



Advantiflex Classic

Premiflex Classic

	HIGH FREQUENCY SEALED AND WELDED SEAMS	INTEGRATED ZIPPER	WATERFALL ZIPPER PROTECTION FLAP	4-WAY STRETCH PATIENT CONTACT SURFACE	COVER WARRANTY	CHLORINE/ SODIUM HYPOCHLORITE RESISTANCE	PREMIUM COMFORT, NON ABRASIVE AND DURABLE UPPER	SUPERIOR DURABILITY, INCREASED SHEAR PROTECTION, ULTIMATE COMFORT UPPER	HEAVY DUTY PVC1100 BASE	HEAVY DUTY PVC1120 BASE	FULL VAPOUR PERMEABILITY	INTEGRATED PRESSURE RELIEF SYSTEM (OPTIONAL)
Advantiflex Classic	✓	✓	✓	✓	2 YEARS	1,000PPM	✓		✓		✓	✓
Premiflex Classic	✓	✓	✓	✓	4 YEARS	10,000PPM		✓		✓	✓	✓



HIGH FREQUENCY SEALED AND WELDED SEAMS

- ✓ Electromagnetically fused seams protecting against risk of penetration of fluid
- ✓ Meets and exceeds infection control requirements
- ✓ Reinforced with additional heavy duty stitching at stress points
- ✓ 100% sealed to prevent ingress of fluid or infection



INTEGRATED ZIPPER

- ✓ Allows easy cover removal
- ✓ Robust and durable
- ✓ High frequency welded attachment for optimal protection against fluid ingress



WATERFALL ZIPPER PROTECTION FLAP

- ✓ Protects against fluid ingress at point of zipper
- ✓ Meets infection control guidelines

ILS INDEPENDENT LIVING SPECIALISTS

1300 008 267

ilsau.com.au



Ph: 1300 731 163

E: sales@fortehealthcare.com.au

W: www.fortehealthcare.com.au