K-Dee II Bed Series







Our most versatile bed offering - it features an attractive, warm design which is lightweight and extremely easy to assemble. Now available in Single or King Single.



Features

The K-Dee II Bed Series is the second generation of the K-Dee Original Bed. It comes fully featured, with powered five-way adjustment, including height, backrest, kneebreak and leg rest. It is compact, lightweight and easy to assemble and disassemble. The Classic Bed features attractive enclosed wooden bed ends and can be fitted with a battery backup. For taller users, an optional length extension can be added to increase the length by 200mm. Both beds are available in either Single or King Single sizes.

- Beech wood grain panels at bed ends on Standard Bed
- Enclosed beech wooden bed ends on Classic Bed
- · Easily assembled and disassembled
- Optional transport bracket for easy storage and transport
- Handset with lockout function
- 4-section pull-apart bed,
- Full power adjustment including Trendelenburg and reverse Trendelenburg
- Dewert actuators

- · Enclosed, individually locking castors
- Dual supports at head end for optional IV pole and self-help pole
- Optional full length wooden or metal rails, or 3/4 collapsible rails
- E1 and Class 3 Fire Retardant
- International standard compliant: BS EN 60601-1:2006 & BS EN 60601-2-252:2010

K-Dee II Standard Bed

Order codes: Single - B09-KDII, King Single - B09-KDII-KS Safe working load: 213kg

Overall length: 2100mm, 2300mm with extension

Overall width: (S) 980mm, (KS) 1160mm

Mattress length: 1980mm, 2180mm with extension

Mattress width: (S) 880mm, (KS) 1050mm Mattress platform height: 335 - 725mm

Backrest adjustment: $0 - 70^{\circ}$ Knee break adjustment: $0 - 23^{\circ}$ Leg rest adjustment: $-27 - 0^{\circ}$

Unit weight: (S) 86kg, (KS) 96kg

Package dimensions: (S) 1100 x 500 x 1200 mm, (KS) 1250

x 550 x 1200mm

Optional accessories: bed extension, IV pole, self-help pole, padded rail covers, transport bracket, full length wooden or

metal rails, 3/4 collapsible rails

K-Dee II Classic Bed

Order codes: Single - B08-KD-CLASSIC, King Single - B08-KD-CLASSIC-KS

Safe working load: 213kg

Overall length: 2150mm, 2350mm with extension

Overall width: (S) 980mm, (KS) 1160mm

Mattress length: 1980mm, 2180mm with extension

Mattress width: (S) 880mm, (KS) 1050mm Mattress platform height: 335 - 725mm

Backrest adjustment: 0 - 70 °

Knee break adjustment: 0 - 23 ° Leg rest adjustment: -27 - 0 °

Unit weight: (S) 90kg, (KS) 100KG

Package dimensions: (S) 1100 x 500 x 1200 mm, (KS) 1250 x

550 x 1200mm

Optional accessories: battery backup, bed extension, IV pole, self-help pole, padded rail covers, transport bracket, full length wooden or metal rails, 3/4 collapsible rails

Notice: all products, product specifications and data are subject to change without notice to improve reliability, function, design or otherwise



sales@cobalthealth.com.au 02 8853 1100 Unit 34, 9 Salisbury Rd, Castle Hill NSW 2154 Australia