

# Donatello Petite Canyon Recliner

Combines modern design aesthetics with an ultra-durable fabric



## 4 Independent motors for a wide range of position combinations and support



### Power Headrest

The individually adjustable head support allows you to lift or lower the position of the headrest even when the backrest is fully reclined providing a higher level of adjustment for your neck and back.



### Power Lumbar

Individually adjustable lumbar support, which can assist in maintaining a proper sitting posture; helping to relieve tension on the spine to prevent pain or discomfort due to poor posture, arthritis or existing back problems.



### Power Backrest and Power Footrest

With a the dual lift system, it allows you to individually adjust the backrest and the footrest for the position that best fits the clinical needs of the client.



## Features

- Dual lift function
- Lay-flat recliner
- Footrest extension
- All close and all lift buttons
- Hand controller
- Ultra durable fabrics
- Lumbar support to accommodate sitting posture

## Easy to use hand control

- Soft grip buttons
- 6 pre-defined positions
- Choose "full recline" or "full stand" with the touch of a button.

General Conditions that the Power Headrest and Lumbar adjustment can assist with:

Kyphosis

Lordosis

Arthritis

Functionality

### Customer Feedback:

I came over to your shop yesterday to try out the latest chair to hit Australia and wow what a chair! There is nothing in your shop the comes close and I've tried them all. Having a fused neck, arthritis and just had two lower back surgeries I am needing support while sitting. I could never get a chair to support my lower back and neck until now because of the electric adjustable lumbar support and electric adjustable head rest. Gone are the days of trying to use pillows to get the support which could never be achieved. Even the way the chair lifts you has improved as you don't feel like you are going to land on your face.

### Therapist Feedback:

I feel this would have a good application for kyphosed clients and offer adjustability not in the typical lift chairs ranges on offer. The contoured seat base is also a good feature, to assist with lower limb and pelvic positioning, acting as a lateral thigh guide to a degree. Although often client with a higher or more complex seating needs may not be the best candidate for the standard lift chairs ranges, these will be welcome additions to provide low to moderate seating support and correction which I think will be beneficial.

### SPECIFICATIONS

Description	Overall Depth	Overall Width	Overall Height	Seat Depth	Seat Width	Seat Height	Reclined Length	Product Weight	Max User Weight	Brand
<b>Donatello Petite Lift Recliner KA55LSC</b>	978mm	851mm	1041mm	457mm	457mm	457mm	1830mm	70kg	158kg	Alivio

Please contact [customerservice@kcare.com.au](mailto:customerservice@kcare.com.au) to find out more information on specific models in each range.  
[www.kcare.com.au](http://www.kcare.com.au)

