


ottobock.

Kidevo adapt

The mobility base for kids.



Max. user weight	75 kg
Max. Payload	100 kg
Seat width range	Size 1: 30 – 40 cm Size 2: 40 – 50 cm Size 3: 50 – 60 cm
Base plate length	Size A: 32 cm Size B: 37 cm Size C: 42 cm
Tilt	-10 to +50°
Drive wheel sizes	20, 22, 24 and 26"
Crash tested (acc. ISO 7176-19)	

ottobock.com.au





Available in three sizes

- The *Adapt* is available in three sizes and allows to adapt the seat width up to 10 cm without any additional parts.



Functional seat tilt

- The *Adapt* has a functional seat tilt ranging from -10° up to $+50^{\circ}$ which allows different positioning.



Reliable anti-tippers

- The anti-tippers support safe driving and maneuvering of the mobility base.



Practical luggage-bin

- Due to the practical luggage-bin, small items can be transported easily.



Angle-adjustable arm support

- The angle-adjustable arm support allows an individual adjustment to the user.



Easy to use brakes

- The *Adapt* can be equipped with an drum brake or combi-drum brake system.



Angle-adjustable push bar

- The push bar is angle-adjustable and provides thereby a comfortable handling.




Variety of adjustment options

- Due to the large variety of adjustment options, optimal driving characteristics can be achieved from active to passive mobility.



Frame colors glossy (RAL-colors)



Frame colors matt (RAL-colors)

