

# ottobock.

## **Kidevo** prime and prime.t

Rigid frame wheelchairs for kids.



Feature	Kidevo prime	Kidevo prime.t
Max. user weight	75 kg	75 kg
Seat width range	26 – 50 cm	28 – 50 cm
Seat depth range	20 – 50 cm	20 – 50 cm
Tilt	–	-10 to + 50°
Drive wheel sizes	22, 24 and 26"	22, 24 and 26"
Crash tested (acc. ISO 7176-19)	✓	✓



### Easy to use brakes

- All *Kidevo* wheelchairs can be chosen with different friendly brakes such as the drum brake. Combinations of brakes are also possible.



### Functional seat tilt

- The *Kidevo prime.t* has a functional and space-saving seat tilt without changing the center of gravity and thereby, maintaining an optimal driving behavior.



### Functional back

- The back of the *Kidevo mini* and *mini.t* is optionally removable, flexible, adjustable in angle and further stabilised.




### Robust foot plate

- For an optimal foot positioning, the *Kidevo* wheelchairs come with a robust and flat one piece foot plate. In addition, there is a high variety of more positioning options to choose from.



### High variety of rear wheels

- Each model of the *Kidevo* provides three different rear wheel sizes in combination with various types of tires and handrims.



### Custom-made seat and back system

- The seat and back system contains of a huge variety of different upholstered custom-made seat and back units.



### Adjustable head supports

- The different types of head supports are adjustable in angle, height and depth.



### High variety of options

- Other options such as belts, therapy trays or grip ring knobs are available. In addition, individual spoke protectors can be chosen for further personalisation.

### Frame colors glossy (RAL-colors)



### Frame colors matt (RAL-colors)



[ottobock.com.au](http://ottobock.com.au)

