



novis molift®
by Etac

Molift Quick Raiser 2

Sit-to-stand hoist with electrically adjustable leg base

Molift Quick Raiser 2 with electrically adjustable leg base is compact, stable and easy to manoeuvre. The unique wheel position and near perfect weight distribution of the base provides a good turning radius and easy manoeuvrability even in narrow spaces.

The inclined straight column imitates the natural pattern of movement and fulfils the hoisting needs of users with some degree of weight bearing capacity. Molift Quick Raiser 2 has a safe working load of 160 kg.

The low base, 10 cm, fits under most beds and chairs while the electrically adjustable legs enables accessibility almost everywhere.



Electrically adjustable legs

A perfect alternative when the hoist is used together with chairs with no accessibility underneath, such as electric wheelchairs and recliners.



Adjustable knee pads

The knee pads are height adjustable and can be turned around for better fitting and to give alternative support for the legs.



Battery indicator

All Molift hoists have an integrated battery indicator that provides audio and visual notification.



Unique wheel positioning

The lightweight construction combined with the unique wheel positioning, provides exceptional maneuverability.



Molift Quick Raiser 2

Item no.

U-shape arm

M04400

Chassis width 610 mm
External width of wheelbase 1200 mm

For complete variant matrix, contact your local distributor.

Total weight

Max user weight (SWL) 160 kg
Chassis weight 36 kg

Material

Aluminum, steel

Dimensions

Height 1280 mm
Leg height 100 mm
Chassis Length 1150 mm

Technical

Battery NiMH 14,4 V - 2,6 Ah
Charge time approx. 3 hrs
Hoisting speed 75 kg 50 mm/sec
Rated performance 40 lifts at 75 kg and 500 mm



Clean / disinfect the product using a standard non-abrasive cleaning agent with a pH between 5 and 9, or with a 70% disinfectant solution.

Proud distributor of Molift products in Australia.

novis

molift®

by Etac

1300 738 885 / novis@novis.com.au

12/12 Mars Rd, Lane Cove West, NSW 2066

Scan QR code for more information

novis.com.au/molift-raiser2

