



**novis molift<sup>®</sup>**  
by Etac

## Molift Smart 150

Compact smartness - light, foldable and transportable without tools

Molift Smart 150 is designed to be easily wheeled and transported to where it is required. It is designed for home care, for travelling and easy storage. The Molift Smart is foldable in one piece without using any tools.

Due to its low weight, small footprint and unique construction, Molift Smart 150 is easy to manoeuvre even in narrow spaces. The hoisting range provides an excellent maximum hoisting height and can also lift from the floor. The 4-point sling bar ensures a comfortable and spacious hoisting position for the user.

Molift Smart 150 has been awarded for Design Excellence by the Norwegian Design Council.



**Foldable without tools**  
Easy handling simplifies everyday life and transport when traveling.



**4-point sling bar**  
Ensures the user a comfortable and spacious hoisting position.



**Easy to manoeuvre**  
Low weight and small footprint facilitates use in narrow spaces.



**Lift from floor**  
Great hoisting interval that allows lifting from the floor.



## Accessories



### Molift Soft Travel Bag

The Molift Soft Travel Bag protects and enables easy transporting of the hoist, battery, charger and sling during travel.



### Molift Travel Suitcase

This hard shell suitcase provides additional protection when transporting by air. It is TSA approved for air travel and ideal for the client with an active life.

#### Product name

Molift Smart 150	M10200
Molift Soft Travel Bag	PHMOLSC
Molift Travel Suitcase	0990200

#### Item no.

#### Technical

Service Software	Included
Battery	NIMH 14.4 V – 2.6 Ah
Battery charge time	approx. 3 hours
Hoisting speed	60 mm/sec
Rated performance	50 hoists/75 kg/500 mm

#### Material

Aluminium, steel

#### Weight

Total Weight	26 kg
Weight when separated	Part 1: 13 kg, Part 2: 13 kg

#### Dimensions

Length	1200 mm
Leg height	110 mm
Dimensions folded	L1160 x W475 x H360 mm
Hoisting height	270–1665 mm
Hoisting interval	1410 mm



Clean / disinfect the product using a standard non-abrasive cleaning agent with a pH between 5 and 9, or with a 70% disinfectant solution. The product can be safely cleaned in a cabinet washer-disinfector at 85°C for 3 minutes.

Proud distributor of Molift products in Australia.

**novis**

**molift**<sup>®</sup>

by Etac

1300 738 885 / novis@novis.com.au

12/12 Mars Rd, Lane Cove West, NSW 2066

Scan QR code for more information

[novis.com.au/molift150](http://novis.com.au/molift150)

