

Oxford[®] Comfort Access Padded Slings




Comfort Access Padded

The Oxford® Comfort Access Padded sling, with head support, has been designed specifically to facilitate the toileting process, allowing easy access to clothing. For use in conjunction with a 4-point positioning cradle, the patented Securi3 clip system ensures no inadvertent detachment of the sling. It is not a general purpose sling and will only suit approximately 25% of clients.


- ✓ Fully padded with integrated head support for additional support
- ✓ Adjustable body belt for added support and security
- ✓ Non-slip back section ensures the sling does not 'ride up'
- ✓ Colour coded binding indicates sling size
- ✓ Positioning handles allow patient repositioning
- ✓ Can be customised to fulfil patient specific needs



Chair to Chair



Chair to Commode/Toilet



Bed to Chair/Chair to Bed
(From seated position only)



Size	Padded with Head Support	SWL
XS	Special Order	227kg
S	LG6394	227kg
M	LG6395	227kg
L	LG6396	227kg
XL	LG6406	227kg





At Enable we source quality, affordable and innovative products globally and have them available for fast delivery. But we do more than just provide products.

We give the knowledge and support to help our customers grow their business and deliver better life care.



ENABLE LIFECARE PTY LTD

T 1300 370 370

enablelifecare.com.au